Neurostimulator

Dolito

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Content

Please read the Operators Manual carefully
before using the device and pay attention
to all references !

	Page
I. Neurostimulator Dolito	
Delivery & Accessories	24
Indication / Contraindication	25
Warnings	26
Advice on Application of Electrodes	27
Technical Data	28
Explanation of Symbols / Maintenance	29
Warranty / Safety Adjustments	30
Real Time Clock Setup	31
II. Display Elements of the Dolito	
Display and Operating Elements	32
Description of the Device	33
III. Application of the Dolito	
Preset Programs	34-35
Setting the Custom Program	36-37
IV. Transcutaneous Electrical Nerve Stimulation (TENS)	
Treatment Modes	38
Programs 01-06 / Abbreviations	39
V. Lock-Button and Statistics	40-41
Placement of Electrodes	42-46
References / Literature	47-48
Waste Disposal / Recycling of Batteries	49
Anatomy	50
Electromagnetic Compatibility	54-56
Copyright	57

Delivery & Accessories

Included in the delivery:

REF: 010E-200: 1 device Dolito, 2 lead wires, 4 self adhesive electrodes (40 x 40 mm), one 9V battery and 1 operators manual

The expected service life of **Dolito** is 10 years.

Accessories:

REF

- 011E-006: 1 set self adhesive electrodes 40x40 (4 pieces)
- 011E-007: 1 set self adhesive electrodes 40x80 (4 pieces)
- 011E-051: 1 set self adhesive electrodes 32mm Ø (4 pieces)
- 011E-008: 1 vaginal probe “IncoProbe V”
- 011E-009: 1 rectal probe “IncoProbe R”
- 011F-012: 2 lead wires
- 011F-013: 2 pieces 9V batteries
- 011F-037: 1 tube electrode jelly for probes
- 011E-208: 1 operators manual Dolito

The battery life is approximately 12 hours or more at 20 mA per channel.

The service life (expiry date) of electrodes, probes and electrode jelly is printed on each package.

Dolito shall only be used with original accessories.

Indications

Dolito was designed for pain-therapy and is easy to use without any assistance:

Nerve Stimulation:

- Increase release of bodily endorphins (PRG 3, 4,)
- Relief or elimination of pain, for instance in the case of:
 - Cervical spine pain (PRG 1, 2, 4)
 - Thoracic spine pain (PRG 1, 2, 4)
 - Lumbar spine pain (PRG 1, 2, 4)
 - Arthritis (PRG 1, 2, 4)
 - Sporting accidents (PRG 1, 2, 4)
 - Circulation problems (PRG 1, 2)
 - Neuralgia (PRG 1, 4)
 - Myalgia (PRG 2, 4)
 - Scar and phantom pain (PRG 2, 3, 4)
 - Fracture pain (PRG 1, 2, 4)
 - Pelvic floor muscle pain (PRG 1, 2, 3, 4)

Contraindications

Dolito may not be used in the following cases:

- Patients with a pace maker as functional disturbances can occur
- The **Dolito** is not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide
- Patients with myocard damages (damage to the cardiac muscle) or arrhythmia (inconsistent heart beat)
- When operating machinery with a potential of danger
- Patients with epilepsy
- When inflammation of the vagina, anus or urinary tract is present (probes)
- Patients with febrile illnesses or infectious diseases

Please pay unconditional attention to the following basic rules:

- Do not stimulate transcerebral (via the head)
- Do not stimulate via the Carotis-Sinus-Nerve (carotid artery)
- Do not stimulate in the area of the Sinus knot (heart)
- Do not stimulate on the Bulbous (eye)
- Do not stimulate in the Larynx-Pharynx-Area (mouth & neck)
- Do not stimulate in places with skin defects or broken skin
- Do not position the electrodes directly on the spine

Warnings

Take special care with the use of the electric stimulation:

- If after an operation the healing process can be interrupted through a muscle contraction
- After acute bleeding or bone fractures with a tendency to bleed
- During menstruation
- In case of desensitized skin with nerve damage

Please note also the following instructions:

If at the same time the connection is made with a surgical high frequency apparatus, this could result in burn marks underneath the electrodes.

If device operation occurs close to short waves or a microwave, instability of the initial values of the stimulator could occur.

The application of the electrodes nearby the thorax can increase the risk of cardiac fibrillation.

In the case of pregnancy, stimulation should only be applied after consultation with a doctor.

If in the case of treatment with the probe, pain or irritation appear in the area of the pelvis, or bleeding in-between menstruation, interrupt the therapy and consult your doctor.

Should skin irritation occur whilst using the **Dolito**, please interrupt the application and inform your doctor - to exclude an allergic reaction.

Patients with metal and/or electronic implants may only stimulate with permission of an experienced doctor.

The **Dolito** may not be operated if the battery compartment is not closed.

Before opening the battery compartment, make sure that there are no cables in the jacks at the head of the unit to eliminate any possibility of current flow.

Never connect the device directly to a battery charger or to any other mains powered equipment.

The operator should not touch the battery contacts & patient simultaneously.

The **Dolito** is protected against debris larger than 12.5 mm in diameter from entering the cabinet and is not water resistant and rated IP20.

Please, keep this device stored out of the reach of children.

Advice on Application of Electrodes

Make sure, that the **Dolito** is switched off (no display visible on the LCD-Display).

Connect the single poled ends of the lead wires to the electrodes.

The output sockets of the electrode cables will be plugged into the channels A or B of the **Dolito**.

Furthermore make sure, that all plugs fit completely into the electrode sockets as well as the device sockets.

Thereafter place the electrodes onto the skin. An adhesive skin gel has been applied to the electrodes, which can be re-used after removal from the skin. In the case of the skin gel losing its adhesion and if you can't achieve more than 8 mA - please replace the electrodes. The same effect, however, also occurs when using a damaged cable. Therefore, if normal treatment is not possible by using new electrodes, the cable maybe faulty.

Electrodes which are damaged due to careless handling may not be used.

Electrodes with damaged cable insulation (e.g. visible copper wire) may not be used.

Electrodes with damaged, missing or stained aluminium coating or carbon layer, also with insufficient adhesive gel may not be used.

Electrodes may not be altered or tampered with (e.g. by scissors or other tools).

Do not pull the electrodes by force on the electrode cable. To remove the electrodes from the skin, please lift the electrode at the rim towards the middle.

Do not use electrodes with ointments or creams (especially which enhance the blood flow or act as local anaesthesia).

Apply the electrodes on clean and healthy skin only.

In the case of bodily hair growth, the hair should be shortened with scissors.

Neither stretch the skin nor the electrodes when applying the electrodes.

Effective current densities of more than 2 mA/cm² require an increased attention of the user with regard to the intensity adjustment.

For standard treatment we recommend the use of our self-adhesive electrodes with an area of 40x40 mm.

If you should use a probe, please clean it carefully before and after use and dry it completely.

Technical Data

Channels:	Two isolated independent channels
Waveform:	Asymmetrical, rectangular. bi-phasic with zero DC current
Intensity:	0 - 80 mA into 500 ohm load 70 mA maximum into 1000 Ohm load 65 mA maximum into 1500 Ohm load (Indication only, actual intensity varies depending on condition of skin electrodes)
Classification:	Internally powered, constant current, continuous operation, Type BF applied part, class IIa (MDD 93/42/EEC), IP20
Frequency:	1 - 150 Hz (Burst 35 - 150 Hz)
Pulse Width:	50 - 350 μ S (Burst 50 - 200 μ S)
Time Settings:	Variable 1 minute to 9 hours 59 minutes (custom program)
Max. Voltage:	180 V without load (70V with 1000 ohm load)
Battery:	9 Volt 6LR61 / 6AM6 recommended

Environmental conditions for storage & transportation:

0-93% atmospheric humidity, -25 to +70 degree centigrade.

Environmental conditions for use: 0-93% atmospheric humidity, +5 to +40 degree centigrade, air pressure 700-1060 hPa.

Dimensions of the device: H 11,7 cm / W 6,6 cm / D 3,3 cm

Weight of the unit: 150 g (incl. battery)

This device is quality controlled and fulfils the rules and regulations set by the German law for medical products (MPG) as well as the guidelines of the European Community 93/42/EEC for medical products. Therefore this device carries the CE-sign „CE 0123“. The serial number is on the rear of the device.

Manufacturer: MTR+ Vertriebs GmbH,
Kamenzer Damm 78, 12249 Berlin - www.mtrplus.com

Explanation of Symbols

	Attention Warning		Manufacturer
	Type BF applied parts		Date of manufacture
	Serial-Number of the product		CE-mark: conforms to essential requirements of the MDD 93/42/EEC
	Keep dry		Not suited for persons with cardiac pacemaker
	LOT number		Read the instructions

Maintenance

Dolito has been designed to be maintenance free. Nevertheless, please adhere to the following advice:

If needed, please clean **Dolito** with a mild detergent, water and a damp cloth. Do not expose **Dolito** to extreme humidity and do not submerge under water. You may clean the cables with a damp cloth. You may also from time to time treat the cables with talcum powder in order to prevent cracks and to further ongoing usage.

A clean cloth, moistened with an alcohol free disinfection liquid can be used to clean the device. Please follow all instructions of the manufacturer of the liquid.

After termination of the therapy, **Dolito** can be re-used by other persons. Please clean it before with a mild disinfectant cleaner.

When changing the battery, please make sure to insert the new batteries within 30 minutes after removal of the empty battery. We recommend that you always have a spare battery with the device.

The **Dolito** does not require calibration.

The **Dolito** is not repairable and contains no user serviceable parts.

No modification to the **Dolito** is allowed.

Warranty

MTR+ Vertriebs GmbH provides a 24 month warranty, commencing from the date of our delivery (invoice date), that in case of failure, the **Dolito** will be serviced or replaced free of charge. All returns must first be authorised by MTR+ Vertriebs GmbH in advance.


Excluded from this warranty are failures of the **Dolito** which result from excessive use, careless or improper handling as well as damage by force.

The warranty loses its lawful meaning if the producer MTR+ Vertriebs GmbH has not given its consent in writing prior to repairworks or servicing being carried out. Without this written consent the manufacturer will not return the device.

The manufacturer reserves the right of technical changes, which may serve to increase treatment success, the functionality of the device or its accessories.

Safety Adjustments

1. The +buttons become inactive after one minute running time of a program to avoid an unintentional increase of intensity. This lock is switched off again by pressing the +button three-times repeatedly. An increasing of intensity is possible again for another minute.
2. If the electrical circuit should be interrupted during operation ($> 8 \text{ mA}$), for example by losing an electrode or by a broken lead wire, the device stops stimulating immediately and the display shows the flashing symbol of an electrode. The effected channel will also flash.

Real Time Clock Setup

Dolito has a “Real Time Clock”, that shows year, month, day, hours and minutes and allows to store detailed statistics of treatment data.

To set time and date, please select the custom program 06 and press the **S**-button for at least 5 seconds. You'll reach the configuration-mode. Now press **S** once again for 10 seconds.

With the +/- button of channel A (left side) you can step forward to the following values

- hour (indicated by: H)
- minute (indicated by: MI)
- day (indicated by: dA)
- month (indicated by: MO)
- year (indicated by: Y)

The respective value is flashing and can be changed with the +/- buttons of channel B (right side). Finally you save all your settings by pressing the PRG-button once again.

Before starting the therapy, please set the clock once if necessary.

If **Dolito** has been without battery for more than 30 minutes, the clock may need to be reset.

Display and Operating Elements

„On/Off“: By pressing the button once, the device will either be switched on or off.

„+/- buttons“: These buttons adjust the intensity for channels A and B. The device has two independently controllable output channels A and B. When pressing the buttons the intensity will either be increased or decreased in steps of 1 mA. 60 seconds after the last change of intensity, a further increase may not be carried out due to an electronic barrier. Pressing the „+“ button three times releases the barrier and the intensity can be altered again. A reduction of the intensity is possible at any time.

Whilst configuring program PC06 one switches between the programmable parameters by pressing the left „+/-“ button (channel A), and through pressing the right „+/-“ button (channel B) the desired value/data for respective parameter can be set. This way, the application time, the stimulation-mode, the frequency and the pulse width can be selected.

„P“: This button selects preset programs 01–05 or custom program PC06.

Press the „P“-button once at any time during a running program to reduce the mA on both channels to zero and pause the program. Pressing „P“ a second time ends the program. Holding „P“ for three seconds while statistics are displayed will exit the statistics screen without deleting them. **Dolito** will return to the home screen of the selected program and the lock will be cancelled.

The „P“ button stores all adjustments made to the custom program mode and returns to the home screen.

„S“: This button is used to setup custom program 06. Press the „S“-button for at least 3 seconds in program PC06 to enter the custom program setting. Pressing „S“ a second time for 10 seconds will enter the real time clock (RTC) setup.

Press the „S“ button once at any time during a running program, to reduce the mA on both channels to zero and pause the program. Pressing „S“ a second time will resume the program and the current will increase over a period of one second back to the previously set intensity level.



Holding the „S“ button for 5 seconds while the statistics are displayed, will display “dEL”.

Holding „S“ a second time for 5 seconds will delete the statistics, set the custom program to default values and show the RTC setup screen. Turn off the unit if you do not want to delete the statistics.

The unit is intended to be supported by hand during normal use. The belt clip is for temporary relief only.

Description of the Device

Preset Programs (01 - 05)

Preparations:

Make sure that **Dolito** is switched off. Cleanse the area of the skin which is to be treated with clean water. Following that, make sure that all plugs are completely connected to the respective electrode sockets as well as device sockets. Place the vaginal-/rectal-probe or the adhesive electrodes onto the area that has been cleansed, according to the advices in pages 25-27 of this manual.


Dolito is operated by a 9 volt block battery. Open the battery compartment at the back by pushing down the lid and insert the battery by hand. Please ensure the correct polarity in accordance with the markings. **Dolito** is operational now.

If the battery voltage is too low, the battery-symbol in the LCD-display will flash. Please replace the battery at this time. Empty batteries are special waste and therefore have to be disposed off at an appropriate collection point or recycled. Never dispose of batteries via domestic or household waste.

Adjustment:


1. Press the button „On/Off” once.

Dolito is now switched on. The program used last, appears on the LCD-display. Program parameters will be displayed as well as the treatment mode and the total treatment time.

2. By pressing the program button „P” you can choose from the programs 01 – 05 (pre-programmed).
3. This is followed by pressing the buttons „+/-” below the channel used. The selected program runs automatically and the LCD displays the intensity as milliamps (mA).

In the upper row you can see the remaining program time and in the bottom line the intensity being applied to each channel. If the button „+” is pressed repeatedly the intensity of stimulation increases. If the button „-” is pressed repeatedly the intensity of stimulation is reduced. For the setting of the intensity it is vital that one feels the electric current but that discomfort never occurs.

4. If you wish to interrupt the program, press the button „S” once and the program will pause. The program will resume for the remaining program time when pressing this button once again. In this case, the intensity increases within one second to the previous used value.
5. By pressing the „P”-button you pause as well, you can continue with new intensity adjustment by pressing the +button. If you press „P” a second time the program is finished.
6. Dolito is switched off when pressing the button „On/Off” or automatically one minute after the program is finished.

Setting the Custom Program (06)

1. Switch on **Dolito** by pressing the „On/Off” button.
2. Choose the program PC06 with the „P” button.
3. Press the „S” button for at least three seconds. **Dolito** changes to the programming mode and the time that indicates program duration, flashes to show that this function can be altered.

In the following please use the left „+“ button to store your choice and to change to the next step. The right „+/-“ buttons are respectively used in order to alter the chosen and flashing parameter.

4. Select with the right „+/-”buttons the desired time for the program.
5. By pressing the left „+”button you will be able to store your choice and change to the next step. Mode now flashes on the LCD and you are able to select the mode that you require by pressing the right „+/-”buttons.

If you have chosen „CONT” or „BURST” you can set the frequency (Hz) and pulse width (μS).

If you have chosen „MOD” you can set the following parameters:

- upper limit of frequency (Hz HI)
- low limit of frequency (Hz LO)
- upper limit of pulse width (μS HI)
- low limit of pulse width (μS LO)

Modulation-time will always be 6 seconds.

Programming of Dense-Disperse-Stimulation according to Prof. Han (**CONT + BURST**) is not possible.

If you have chosen “WORK/REST”- rhythm you can set the following parameters.:

- frequency (Hz)
- pulse width (μS)
- work time (**WORK**)
- rest time (**REST** = recovery time)

6. The programming is completed by pressing the „P” button and the program is stored.

Selectable frequency and pulse width settings:

Continuous-Mode: 1 to 150 Hz, 50 to 350 μ S

Burst-Mode: 35 to 150 Hz, 50 to 200 μ S

Modulated-Mode: 1 to 150 Hz, 50 to 350 μ S

Work/Rest-Mode: 1 to 150 Hz, 50 to 350 μ S

The settings can be made as follows.

- Frequency: 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 15 Hz, and then increasing in 5 Hz steps until max. 150 Hz
- Pulse width: 50 to 350 μ S in steps of 10 microseconds
- Work time: 2-99 seconds in steps of 1 second
- Rest time: 2-99 seconds in steps of 1 second
- Total treatment time: 1-599 minutes in steps of 1 minute

TENS (Transcutaneous Electrical Nerve Stimulation)

TENS as a scientifically recognised method for pain relief therapy is a suitable aid for pain for many doctors. The stimulating electrical current impulses being produced by the **Dolito** nerve stimulator are transferred from the electrodes to the skins nerve fibres. The effect of pain relief begins at the start of the application and is normally effective for an amount of time after the stimulation period. This is achieved due to two different mechanisms:

Gate-Control-Effect:

The pain is transferred over the nerve fibres of the body to the brain which perceives the pain. Through stimulus of the nerves the **Dolito** nerve stimulator prevents further transmission of the pain signals to the brain. The pain sensation is therefore eased or inhibited.

Beta-endorphins-release:

TENS-therapy stimulates the synthesis and release of beta-endorphins. Beta-endorphins as bodily synthesised peptides, which have a similar effect to morphine for pain relief, but without the side effects that morphine carries (e.g. breathing depression).

TENS is available in three different **treatment modes**.

- 1. Continuous Tens (CONT):** Electrical impulses are exchanged between the electrodes continuously back and forth. They block the transmitting nerve bands. It is possible that the body gets used rather quickly to this traditional treatment. The effects can subside with prolonged treatment.
- 1a** A special kind of continuous mode is named **Dense-Disperse-Stimulation (CONT + BURST)** according to Prof. Han. The permanent alternation of frequency between 2 Hz and 100 Hz within a few seconds (2-4 secs.) raises the output of different analgesic substances.
- 2. Modulated Tens (MOD):** By changing the parameters (frequency, pulse width and/or intensity) regularly and automatically, the sensation of an increasing and decreasing current is produced. Therefore the body will adapt distinctly delayed to the electrical impulses.
- 3. Burst (BST):** 9 impulses with a set frequency and pulse range, twice per second are emitted. This artificially produced frequency of 2 Hz is able to increase the production of the bodies own endorphins and therefore relieve pain in a more natural way.

Programs

Program	01	02	03	04	05	06
Mode	MOD	CONT	BURST	HAN	W/R	selectable
Frequency	70-110 Hz	110 Hz	80 Hz	100/2 Hz	50 Hz	1-150 Hz
Pulse Width	200-100 μ S	175 μ S	150 μ S	200/350 μ S	300 μ S	50-350 μ S
Running Time	30 min.	30 min.	30 min.	30 min.	20 min.	1-599 min.
Ramp up					1 sec.	1 sec.
Ramp down					1 sec.	1 sec.
Work Period					10 sec.	2-99 sec.
Rest Period					20 sec.	2-99 sec.

Abbreviations:


- AdJ: Adjustment of intensity
AL: Alternate function of the channels A and B
BURST: The burst mode delivers 2 blocks of 9 impulses each per second
CONT: In the continuous mode a permanent unchanged current is used
dA: Day
dEL: Delete
dL: Delay between channels A and B
HI: Upper value (high)
Hz: Hertz is the measurement unit for the frequency of impulses per second.
LO: Lower value
mA: Milli-amp is the measurement unit for the electrical current intensity.
MOD: In the modulated mode the frequency and/or the pulse width and/or the intensity changes in a predefined rhythm automatically.
Rd: Ramp down
RU: Ramp up
SEC: Seconds
SY: Synchronous function of the channels A and B
TI: Time
TM: Time of modulation-interval
W/R: Work/Rest mode
 μ S: Microsecond is the measurement unit for the length of every single impulse.

LOCK-Button and Statistics


The **Dolito** offers the ability to lock one program in place so that other programs or settings cannot be selected or changed. This allows the doctor or therapist to ensure that the patient uses only the program that has been selected for him. The therapist can read information about the compliance of the patient on the display at any time he wants. This way it is possible to evaluate, whether the device was used according to the instructions of the doctor (quantity, intensity, etc.)

For this purpose the LOCK-button in the battery-compartment has to be pressed with a small pin once (the metal-pins of the electrode-cables are well suited). The lock symbol appears on the display and the program cannot be changed by pressing the „P“⁺-button anymore. Pressing the LOCK-button a second time allows reading the following statistics for the treatment in the display:

If the program was not locked, it is necessary to press the LOCK button twice to read the memory. The statistics are always recorded regardless of whether the unit is locked or unlocked. The statistics should be reset after 60 days of use or if the device is used by another person.


To get more detailed information press the + button of the left channel (A). You will see the data of the first session of the first day.

The next sessions for the selected day are displayed - one after the other - by pressing the „S”- button. The next day can be obtained by pressing the + button of channel A again. Hold A+ or A- for three seconds to return to the global statistics screen.

If you require information about time and date of one single session please press any of the right +/-buttons (channel B). With these buttons you can scroll between session data, date and time.

The data stored in memory is retained after viewing, so that therapy-control is possible anytime. A maximum of 60 days of data can be stored with 5 sessions per day. When 60 days of data are reached, consecutive days will overwrite the earliest days. The day number will continue increasing (61, 62, ... up to max. 180) and will be displayed when the unit is unlocked until the **Dolito** is reset.

If you press the P button for three seconds, you return to the home screen of the selected program. The device will not be locked. Press the lock button once to lock the device if it is required.

If you do not press „P” for 3 seconds, you can turn off the unit to leave it in lock mode when it is restarted.

The statistics should be reset before a new patient uses the device. Press LOCK to enter the statistics display and press „S” twice for 5 seconds each time. All settings in program PC06 will be reset as well. The **Dolito** is set to its default state and all statistics are deleted.

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In most cases it is very easy to find the suitable points for the placing of the electrodes. It is recommended to seek advice on electrode placement from your doctor or physiotherapist. We suggest placing the electrodes on either side of the pain if no advice is available. The distance between the electrodes should not be more then 20 cm. The polarity of the electrodes has no significance in TENS therapy.


Examples of possible placement-points are shown on the following pages.

While stimulating, please keep this device in reach, to be able to finish the treatment at once, if current becomes uncomfortable.


Please ensure you follow the warnings and electrode placement advice on pages 25-27 of this manual.

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















Menstruationsschmerzen
Menstrual Pain

Literatur Literature

1. **Benton, L. / Baker, L.L. / Bowman, B.R. / Waters, R.L.:** Funktionelle Elektrostimulation. (1983) Steinkopff Verlag Darmstadt
2. **Birklein, F.:** Neuropathische Schmerzen – Mechanismus und Therapie, Fortschr Neurol Psychiat; 2002, 70: 88-94
3. **Bjordal, J.M. / Johnson, M.I. / Ljunggreen, A.E.:** Transcutaneous electrical nerve stimulation (TENS) can reduce postoperative analgesic consumption. A meta-analysis with assessment of optimal treatment parameters for postoperative pain. European Journal of Pain 7 (2003), 181-188
4. **Bossert, F.-P. / Vogedes, K.:** Elektrotherapie, Licht- und Strahlentherapie, Urban & Fischer Verlag, Jena/München 2003
5. **Chen, X.H. / Han, J.S.:** Analgesia induced by electroacupuncture of different frequencies mediated by different types of opioid receptors: another cross-tolerance study; Behavioural Brain Research, 47 (1992) 143-149
6. **Fisher, W.W. / Bowman, L.G. / Thompson, R.H. / Contrucci, S.A. / Burd, L. / Alon, G.:** Reduction in self-injury produced by transcutaneous electrical nerve stimulation, Journal of Applied Behavior Analysis 31 (1998), 493-496
7. **Gillert, O. / Rulffs, W. / Boegelein, K.:** Elektrotherapie 3. Auflage (1995) Pflaum-Verlag, München
8. **Gupta, A.K. / Gupta, S. / Meena, D.S. / Sharma, U.:** Post-Tonsillectomy Pain: Different Modes of Pain Relief, Indian Journal of Otolaryngology and Head and Neck Surgery. Vol. 54, No. 2, (2002)
9. **Halbert, J. / Crotty, M. / Cameron, I.D.:** Evidence for the Optimal Management of Acute and Chronic Phantom Pain: A Systematic Review, The Clinical Journal of Pain 18 (2002), 84-92
10. **Hamza Mohamed A. / White, Paul F. / Hesham, E. Ahmed / El-Sayed, A. Ghoname:** Effect of the Frequency of Transcutaneous Electrical Nerve Stimulation on the Postoperative Opioid Analgesic Requirement and Recovery Profile, Anesthesiology (1999); 91: 1232-8
11. **Han, J.S.:** Acupuncture: neuropeptide release produced by electrical stimulation of different frequencies; Trends in Neurosciences, Vol.26, No.1, January 2003
12. **Han, J.S. / Chen, X.H. / Sun, S.L. / Xu, X.J. / Yuan, Y. / Yan, S.C. / Hao, J.X. / Terenius, L.:** Effect of low- and high-frequency TENS on Met-enkephalin-Arg-Phe and dynorphin A immunoreactivity in human lumbar CSF; Pain, 47 (1991) 295-298
13. **Jenrich, W.:** Grundlagen der Elektrotherapie, Urban & Fischer Verlag, Jena/München 2000
14. **Johnson, M.:** Transcutaneous Electrical Nerve Stimulation (TENS). In Kitchen, S. ed.: Electrotherapy – evidence based practice; Edinburgh; Churchill Livingstone, 2001: 259-286
15. **Johnson, M.I.:** Transcutaneous Electrical Nerve Stimulation (TENS) and TENS-like devices: do they provide pain relief? Pain Reviews 8 (2001); 7-44

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In accordance with European regulations [1], used electrical and electronic equipment may not be discarded in house hold or unsorted municipal waste. Use separate collection facilities. Contact your local government or council for information regarding the collection systems available. If electrical appliances are disposed of in landfills or dumps, hazardous substances can leak into the ground water and get into the food chain damaging the health and well being of people and animals.

[1] European Community directive 2002/96/EC on waste electrical and electronic equipment.

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
The battery should be replaced when it is drained or stops functioning. The battery supplied with the Dolito is recyclable. Remove the battery from the unit and follow your local recycling guidelines.

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Electromagnetic Compatibility

Guidance and manufacturer's declaration – electromagnetic emissions		
The Dolito unit is intended for use in the electromagnetic environment specified below. The customer or the user of the Dolito should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Dolito uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Dolito is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not Applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable	

Guidance and manufacturer's declaration – electromagnetic immunity			
The Dolito unit is intended for use in the electromagnetic environment specified below. The customer or the user of the Dolito should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/-6 kV contact +/-8 kV air	+/-6 kV contact +/-8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Electromagnetic Compatibility

Guidance and manufacturer's declaration – electromagnetic immunity			
The Dolito unit is intended for use in the electromagnetic environment specified below. The customer or the user of the Dolito should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 V_{rms} 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2,5 GHz</p>	<p>3 V</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Dolito unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p> <p>$d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ 80MHz to 800MHz $d=2.3 \sqrt{P}$ 800MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: right;">  </div>
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>A) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Dolito unit is used exceeds the applicable RF compliance level above, the Dolito unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Dolito unit.</p> <p>B) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.</p>			

Electromagnetic Compatibility

Recommended separation distances between portable and mobile RF communications equipment and the Dolito unit.			
<p>The Dolito unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Dolito unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Dolito unit as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
	$d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$	$d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$	$d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

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Reserved for changes.

Printing or duplication of this manual - also in an abridged version - is forbidden and can only be allowed with a written approval of “MTR+ Vertriebs GmbH”.

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