Stimulator für paretische Muskulatur
TENS & EMS
Stimulator for paretic muscles
TENS & EMS

PeroBravo

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Content

Please read the Operators Manual carefully before using the device and pay attention to all references !

Page

I. Stimulator for Nerves and Muscles (incl. Paresis) PeroBravo

Delivery & Accessories / Explanation of Symbols	34
Technical Data	35
Introduction - What is PeroBravo	36
Indications / Contraindications	37
Contraindications / Safety Adjustments	38
Warnings	39
Advice on Application of Electrodes	40-41
Accessories	42
Maintenance / Warranty	43

II. Operating- and Display-Elements of PeroBravo

Operating Elements / Connections	44
Battery Compartment with LOCK-button	45
The Display	46
The Info-line	47
Real Time Clock Setup / Language Selection	48

III. The Programs of PeroBravo

Single Pulses	49
Biphasic Current	50-51
Monophasic Current	51
Diagnostic	52
Walking Training	53
Walking Training - Custom Program	54
Single Pulses - Custom Programs	55
Biphasic Current / Monophasic Current - Custom Programs	56-57

Content

	Page
IV. Foot-Switch	58
V. Electrode Placement	
Electrode Placement Paresis	59
Electrode Placement EMS	61-67
Electrode Placement TENS	68-72
Anatomy	73
VI. Control of Therapy	60
VII. Batteries and Battery-Charger	74
Waste Disposal / Recycling of Batteries	75
References / Literature	76-77
Electromagnetic Compatibility	81-83
Copyright	83
SD-curve	84-85

Delivery & Accessories

Included in the package:

REF: 010E-112: 1 unit **PeroBravo**,

4 self-adhesive electrodes (40 x 40 mm), 1 cable to connect the electrodes, 1 foot-switch, 2 rubber electrodes (45 x 50 mm), 2 sponge-bags, 2 straps to fix the rubber electrodes, 4 rechargeable batteries type AA, 1 battery-charger, 1 carrying case and 1 operators manual

Accessories:

REF

011E-006: 1 plastic bag self adhesive electrodes 40x40 (4 pieces)

011E-007: 1 plastic bag self adhesive electrodes 40x80 (4 pieces)

011E-113: 1 foot-switch

011E-114: 2 rubber electrodes 45 x 50 mm

011E-115: 2 sponge bags

011F-116: 2 velcro straps

011F-117: 4 rechargeable batteries type AA

011F-118: 1 battery charger

011E-119: 1 carrying case

011E-120: 1 cable to connect the electrodes

011E-121: 1 operators manual

Explanation of Symbols

Attention
Warning

Type BF applied parts

Serial-Number
of the product

Keep dry

Batch code

Indoor use only

Manufacturer

Date of
manufacture

CE-mark: conforms to
essential requirements of
the MDD 93/42/EEC

Not suited for persons
with cardiac pacemaker

Protection class II

Read the instructions

Technical Data

- 1-channel stimulator with 20 programs.
- Intensity 0 - max. 60 mA into 500 ohm or 1000 ohm load, 50 mA maximum into 1500 ohm load (Indication only, actual intensity varies depending on condition of the electrodes) With monophasic current max. intensity is 40 mA.
- Power supply: 4 rechargeable (NiMH-)batteries 1,5 V Type AA
- Voltage max.:130 V without load (70 V into 1000 ohm load)
- Dimensions: L 141,6 mm x W 74 mm x H 40,4 mm
- Weight: 275 g (incl. batteries)

Single pulses:

- Waveform: monophasic triangular or trapezoidal with continuously changing polarity
- Pulse width triangle: 1 ms - 1000 ms / Pulse width trapeze: 5 ms - 1000 ms
- PRG 12: Pulse width monophasic rectangle: 1 ms - 1000 ms

Enveloped pulses:

- Waveform: biphasic, symmetrical, rectangular or monophasic positive, rectangular
- Frequency: 1 - 120 Hz
- Pulse width: 50 μ S - 35 ms (positive part)

The **PeroBravo** features a stimulation circuit failure detection system incorporated into the device and its software. Four pulses with a pulse width and voltage not exceeding 1 μ S and 6 volts are triggered before each main pulse to check the integrity of the circuit. These pulses cannot be felt and do not effect the efficacy of the device. If a failure is detected during this test, STIM.FAIL will display on the LCD. If STIM.FAIL is repeatedly displayed on the LCD then the unit should be returned to the manufacturer for testing. This is an added safety protection system to ensure each pulse generated by the device is safe for the user.

Classification: Internally powered, constant current, continuous operation, Type BF applied part, class IIa (MDD 93/42/EEC), IP 20

The expected service life of **PeroBravo** is 10 years. Service life of accessories is one year. The storage time of the batteries is approximately 2 years and the charging life is about 500 cycles. The service life (expiry date) of self-adhesive electrodes is printed on each package.

Environmental conditions for storage & transportation: 0-93% atmospheric humidity, -25 to +70 degree centigrade. Environmental conditions for use: 0-93% atmospheric humidity, +5 to +40 degree centigrade, air pressure 700-1060 hPa.

This device is quality controlled and fulfils the rules and regulations set by the German law for medical products (MPG) as well as the guidelines of the European Community 93/42/EEC for medical products. Therefore this device carries the CE-sign „CE 0123“. The serial number is in the battery compartment.

Manufacturer: MTR+ Vertriebs GmbH,
Kamenzer Damm 78, 12249 Berlin - www.mtrplus.com

Introduction - What is PeroBravo

PeroBravo is a versatile stimulator with 20 different programs, which can be used by the patient without any help. The unit is intended to be supported by hand during normal use. The belt clip is for temporary relief only. The device uses monophasic or biphasic currents both for the treatment of pain as well as for muscle stimulation. One focus is the treatment of paresis, especially peroneal palsy.

If any patient suffers from paresis the nerves leading to a muscle or group of muscles are disturbed for various reasons. The voluntary activation of the muscle is no longer possible and the muscle is inactive. Due to this inactivity, the number of active muscle fibers will reduce and the muscle mass will decrease within a few weeks (atrophy). Electrotherapy will stop this atrophy or at least slow down, so that in case of successful regeneration of the nerves the muscle can still be functional.

The treatment of paresis by means of electric stimulation contains several Types of treatment: The initial treatment usually consists of long single pulses in a triangular or trapezoidal shape, as the paretic muscle does not respond to the usual short rectangular pulses of normal muscle stimulation. In the regeneration phase of the nerve muscular atrophy shall be prevented and regeneration shall be accelerated by repetitions of short single twitches. This is usually done with pulses using a long ascent. The length of ascent ensures a reaction of the paretic muscle and prevents the surrounding healthy muscles from contraction as well. For less serious conditions, or with advancing therapeutic success, the pulses may be shorter and/or can be modified from the triangular shape towards a trapeze - but only so far that the partial denervated muscle still responds to the stimulus. The optimal treatment parameters can be determined individually by creation of an strength/duration-curve (SDC) for each patient.

One of the most common kind of palsy is peroneus-paralysis or foot-dorsiflexion. With this disease the patient is no longer able to raise the toe of the affected leg. If the muscle is ready to respond on rectangular pulses, it is possible to start the second phase of treatment with normal muscle training by short biphasic rectangular pulses in a certain frequency and with work- and rest-periods. At this stage, training of walking can start as a third step. It is supported by a foot-switch that triggers current for the swing-phase of the gait.

With the **PeroBravo**, a therapy device is available now that combines all treatment sections in one device. We tried to create a new stimulator that helps to reduce cost for health service without limitations regarding the quality of treatment.

Indications

PeroBravo is mainly applicable for the treatment for all forms of flaccid paralysis. This applies both for the initial treatment of severe paresis (PRG 1-6, 17, 18), as well as for treatment in an already improved stage (PRG 7-11, 19, 20). A special program allows the creation of an I/t-curve (PRG 12) to find the optimal treatment parameters for each patient individually.

In addition, the included footswitch allows gait-synchronous training of walking for patients with peroneal palsy using biphasic rectangular pulses in a certain frequency (PRG 13-16).

The device is not only the first choice for partial denervated muscles, but can also be used for muscle rehabilitation with intact nervous control (PRG 10, 19, 20) or for pain management (PRG 19, 20).

Since the programs 19 and 20 can be defined according to individual requirements **PeroBravo** can be used comfortably for the following treatments:

Muscle Stimulation:

- Avoidance and delay of muscular spasms
- Relaxation of muscular spasms
- Avoidance and delay of muscular amyotrophy
- Recovery of the muscle
- Conservation and extension of mobility
- Extension of the capillary bed
- Speeding up the process of removal of toxic substances from metabolism
- Strengthening of the muscle
- Activation of former inactive muscular fibres
- Increase resistance of the muscle against fatigue
- Avoidance of muscular imbalances

Nerve Stimulation:

- Increase release of bodily endorphins
- Relief or elimination of pain, for instance in the case of:
 - Cervical spine pain
 - Thoracic spine pain
 - Lumbar spine pain
 - Arthritis
 - Sporting accidents
 - Circulation problems
 - Neuralgia
 - Myalgia
 - Scar and phantom pain
 - Fracture pain
 - Pelvic floor muscle pain

Contraindications

PeroBravo may under no circumstances be used in the following cases:

- Patients with electronic implants like pumps or a pace maker because functional disturbances can occur
- The **PeroBravo** is not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.
- Patients with myocard damages (damage to the cardiac muscle) or arrhythmia (inconsistent heart beat)
- When operating machinery with a potential of danger
- Patients with epilepsy
- Patients with febrile illnesses or infectious diseases

Please pay unconditional attention to the following fundamental rules:

- Do not stimulate transcerebral (via the head).
- Do not stimulate via the Carotis-Sinus-Nerve (carotid artery).
- Do not stimulate in the area of the Sinus knot (heart).
- Do not stimulate on the Bulbous (eye).
- Do not stimulate in the Larynx-Pharynx-Area (larynx /pharynx).
- Do not stimulate in places with skin defects.
- Do not position the electrodes directly on the spine.

Take special care with the use of the electric stimulation:

- if after an operation the healing process can be interrupted through a muscle contraction.
- after acute bleeding or bone fractures with a tendency to bleed.
- in case of desensitized skin with nerve damage.
- if stimulation is located close to osteosynthesis material (metal in the body).
- if patients cannot use the device properly, due to mental or physical limitations.

Safety Adjustments

1. The +buttons become inactive after one minute running time of a program to avoid an unintentional increase of intensity . This lock is switched off again by pressing the +button three-times repeatedly. An increasing of intensity is possible again for another minute.
2. If the electrical circuit should be interrupted during operation, for example by losing an electrode or by a broken lead wire, the device stops stimulating immediately and the display shows the flashing symbol of an electrode. The effected channel will also flash.
3. If any phase comes to its end, the intensity falls to 8 mA automatically and must be increased by the user to the desired intensity once again.
4. Between two single pulses, intensity cannot be increased more than 3 mA.
5. If the area under the trapezoidal single pulses increases by shortening the ascent, the intensity is automatically reduced proportionally.
6. Maximum intensity in each program is limited according to frequency and pulse width. The limit in monophasic programs is max. 40 mA.

Warnings

PeroBravo is a type BF equipment and may not be immersed in liquids.

PeroBravo may not be used with a vaginal or rectal probe.

If at the same time the connection is made with a surgical high frequency apparatus, this could result in burn marks underneath the electrodes.

If device operation occurs close to short waves or a microwave, instability of the initial values of the stimulator could occur.

The application of the electrodes nearby the thorax can increase the risk of cardiac fibrillation. In the case of pregnancy, stimulation should only be applied after consultation with a general practitioner.

Should skin irritation occur whilst using the **PeroBravo**, please interrupt the application and inform your doctor - to exclude an allergic reaction.

Patients with metal and/or electronic implants may only stimulate with permission of an experienced doctor.

In the area of tumors stimulation is allowed only by instruction of the specialized physician.

If high frequency is combined with large pulse width, please use only small intensity (see safety limits on page 57).

The **PeroBravo** may not be operated if the battery compartment is not closed. Before opening the battery compartment, make sure that there are no cables in the jacks at the head of the unit to eliminate any possibility of current flow.

If using rechargeable batteries be sure to use a CE approved battery charger.

Never connect the device directly to a battery charger or to any other mains powered equipment.

The operator should not touch the battery contacts & patient simultaneously.

The **PeroBravo** is protected against debris larger than 12.5 mm in diameter from entering the cabinet, is not water resistant and rated IP20.

If you have swallowed any part of the device or any of its accessories, please contact a doctor immediately.

Electrodes can contain amounts of nickel. Allergy sufferers should check this before using.

Please, keep this device and its accessories stored out of the reach of children.

Advice on Application of Electrodes

Make sure, that the **PeroBravo** is switched off (no display visible on the LCD-Display).
Connect the single poled ends of the lead wires to the electrodes plugs.
The other end of the electrode cable will be plugged into the output jack A of the **PeroBravo**.
Make sure, that all plugs fit completely into the electrode sockets as well as the device socket.
Thereafter place the electrodes onto the skin.

Self-adhesive electrodes:

The self-adhesive electrodes are intended for gait training in the programs 13-16 and for use with biphasic rectangular pulses that are shorter than 500 μ S. An adhesive skin gel has been applied to the electrodes, which can be re-used after removal from the skin. In the case of the skin gel loosing its adhesion and if you can't achieve more than 1 mA - please replace the electrodes. The same effect, however, also occurs when using a damaged cable. Therefore, if normal treatment is not possible by using new electrodes, the cable maybe faulty.

Electrodes which are damaged due to careless handling may not be used.

Electrodes with damaged cable insulation (e.g. visible copper wire) may not be used.

Electrodes with damaged, missing or stained aluminium coating or carbon layer, also with insufficient adhesive gel may not be used.

Don't stretch the skin or the electrodes when applying the electrodes.

Do not pull the electrodes by force on the electrode cable. To remove the electrodes from the skin, please lift the electrode at the rim towards the middle.

After use, please put back the self-adhesive electrodes on the foil and keep it in the sealed plastic bag in a cool place.

These electrodes are intended for use on one single patient only.

Rubber electrodes:

Rubber electrodes must be used in all programs with monophasic pulses or pulses with a pulse width longer than 500 μ S.

Please moisten the sponge bags with water thoroughly. Then make sure the rubber electrodes are fully inserted into the wet sponge pockets.

Now fix the sponge-bags with electrodes to the desired locations with the help of the velcro straps.

Make sure that the contacts of the connection cables are fully inserted into the sockets of the electrodes before turning on the device. If the electrode contact point has not been met exactly, you can easily move the sponge bags under the straps to the right position - to optimize the stimulation.

After each treatment, please take the rubber electrodes out of the sponge bags, disconnect the cables and rinse the electrodes as well as the sponge bags with clean water. Then the rubber electrodes should be dried carefully.

Advice on Application of Electrodes

In general:

Electrodes may not be altered or tampered with (e.g. by scissors or other tools).

Do not pull the electrodes by force on the electrode cable.

Do not use electrodes with ointments or creams (especially which enhance the blood flow or act as local anaesthesia).

Apply the electrodes on clean and healthy skin only.

In the case of bodily hair growth, the hair should be shortened with scissors.

Effective current densities of more than 2 mA/cm² require an increased attention of the user with regard to the intensity adjustment. Therefore we either recommend the use of our self-adhesive electrodes with a minimum area of 40x40 mm or the use of rubber electrodes with wet sponge bags. Maximum recommended intensity is 53 mA for 40x40 mm electrodes. Electrodes with 45x50 mm or 40x80 mm may be used up to 60 mA.

All custom electrodes with the above minimum size which comply with the requirements set out in the Council Directive 93/42/EEC relating to the Medical Device Directive can be used, if the connections match the cable of **PeroBravo**.

Pulses over 500 µS and monophasic pulses can increase the risk of skin irritation. Using rubber electrodes with wet sponge bags reduces this irritation.

Always use rubber-electrodes with wet sponge bags in all programs which work with monophasic pulses (1-6, 11, 12, 17, 18, 20) - even if these pulses have changing direction.

All materials that have been brought to the skin by the sponge-bags, electrodes or velcro straps can be removed completely with water and soap.

Accessories

Sponge bags:

The sponge bags are made from commercial viscose and consist of a mixture of cellulose and cotton. The material is kind to the skin and meets the requirements of DIN EN 13432 (compostable material).

Before first use, please rinse with warm water once.

The material provides a high water absorption. Please make sure that the treatment is done with a wet sponge only. If the sponges are slightly smaller and stiffer after drying, this does not effect the quality. Once they are soaked for therapy again, they get back their original flexibility.

If necessary, the sponge bags can be cleaned with clear, warm water up to 95°. Please avoid the use of detergents or disinfectants, ironing, bleaching, chemical cleaning and tumble drying.

Please airdry the sponge bags after use completely, don't put them back into the plastic bag when damp!

These sponge bags are intended for single patient use only.

Elastic velcro straps:

The elastic velcro straps are used for rapid attachment of rubber-electrodes and sponge-bags. They allow the sponge bags to move under the straps to an optimum position and provide comfort. Please do not over tighten the straps as they may interfere with blood circulation.

Under normal use the velcro straps are durable and can be washed several times in warm water with a temperature of up to 60 ° and a standard household laundry detergent. Please close the velcro fasteners before washing and avoid ironing, bleaching, chemical cleaning and tumble drying.

These straps are intended for single patient use only.

Maintenance

PeroBravo has been designed to be maintenance free. Nevertheless, please adhere to the following advice:

If needed, please clean **PeroBravo** with a mild detergent, water and a damp cloth. Do not expose **PeroBravo** to extreme humidity and do not submerge under water. You may clean the cables with a damp cloth. You may also from time to time treat the cables with talcum powder in order to prevent cracks and to further on going usage.

After termination of the therapy, **PeroBravo** can be re-used by other patients. Please clean it before with a mild disinfectant cleaner.

Please subject the device to technical safety control in regular intervals in accordance with the local regulations. MTR+ recommends servicing the unit every two years. The **PeroBravo** does not require calibration. The **PeroBravo** is not repairable and contains no user serveable parts. No modification to the **PeroBravo** is allowed.

Remove the batteries from the **PeroBravo** if it is not used for extended periods of time.

Warranty

MTR+ Vertriebs GmbH provides a 24 month warranty, commencing from the date of our delivery (invoice date), in case of failure, the **PeroBravo** will be serviced or replaced free of charge. All returns must first be authorised by MTR+ Vertriebs GmbH in advance.

Excluded from this warranty are failures of the **PeroBravo** which result from excessive use, careless or improper handling as well as damage by force.


The warranty loses its lawful meaning if the producer MTR+ Vertriebs GmbH has not given its consent in writing prior to repairworks or servicing being carried out. Without this written consent the manufacturer will not return the device.

The manufacturer reserves the right of technical changes, which may serve to increase treatment success, the functionality of the device or its accessories.


Operating Elements

The SET button allows access to the program parameters that are individually selectable when it is held down for at least 5 seconds. A single press during a program pauses the program, it is resumed on the second press with the previously set intensity. This function is disabled for the gait training (PRG 13-16). If “manual trigger” was selected in the programs 1-6 or 17-18 (see page 49, 55), the pulse can be triggered manually here.

The PRG-button is used for program-selection. During a running program, pressing it will pause the program. A second press will finish the program. In addition custom parameters are stored with this button.

Connections


The programs 13-16 are available only when the foot-switch is connected. All other programs are not accessible while the footswitch is connected.

Battery Compartment with LOCK-button


PeroBravo allows the doctor or therapist to lock one program for the patient to use. In this way the medical practitioner can be quite sure that the patient only uses the device with the correct program. The LOCK-button can be pressed with a pin of the cable or any other suited thin object. Pressing will fix the selected program so that different programs cannot be used anymore.

UNLOCK-button removes the program lock and shows global statistics

Battery orientation


At the front of the battery compartment you will find the serial number and production date next to the barcode.

The Display

The screen of **PeroBravo** consists of two different parts. The first part - above and below - includes fixed signs and symbols. All important informations for operation and therapy-monitoring are indicated here by numbers and letters. The second - shaded - area in the middle provides an optical impression of the waveform that is currently used in the particular program and other useful information.

The Info-Line

The info-line in the upper part of the screen shows the remaining treatment-time if programs 7-11 or 19, 20 are selected. In programs with single pulses the percentage of ascent, top level and descent of the pulse is displayed.

Outside of a treatment session the info-line shows the following informations:

Setup real time clock: SET HR - SET MIN - SET DAY - SET MTH - SET YR

Treatment modes:

CON (CONSTANT) = Continuous mode, no automatic changes in pulse width, frequency or intensity

MOD (MODULATE.) = Modulated mode, periodic automatic changes in pulse width, frequency and/or intensity

BURST = Bundling of each 9 pulses to pulse groups twice per second

W/R (WORK.REST) = Work/Rest-rhythm of contraction and relaxation in muscle training

RIC = Rectangle Impulse Characteristic for creation of a strength/duration-curve

TIC = Triangle Impulse Characteristic for creation of a strength-duration-curve

Shape of the pulses: TRIANGLE - TRAPEZE

Notes on certain parameters (for programming):

PULSE = Adjust pulse width

PAUSE = Define pause-time between single pulses / Interruption of program

REPS = Define number of repetitions

TRIGGER.A / TRIGGER.M = Select automatic or manual trigger of pulses

TIME = Define length of the phase

NO LIMIT = Duration of phase or program unlimited

FREQ. = Define frequency in programs 7-11, 13-16, 19, 20

HI / LOW = Define upper or lower value

INTEN.MOD. = Intensity modulation selectable from 100% (no modulation) up to 50%

TIME MOD = Define the time of one modulation interval

WORK = Select duration of the muscle contraction

REST = Select duration of the relaxation between the contractions

RAMP UP = Duration of the ramp up to 100% of the selected intensity

RAMP DN = Duration of the ramp down at the end of each contraction

CONNECT = Prompts for connection of the foot-switch when using programs 13-16

SW HEEL = Mode used if foot-switch is placed under the heel (standard)

SW TOE = Mode used if foot-switch is placed under the healthy foot (exceptional option)

START.DLY = Sets a delay time for the start of stimulation to optimize the foot in gait

STOP.DLY. = Sets a delay time for the end of stimulation to optimize the foot in gait

PRG END = End of program

CABLE = Current circuit is interrupted by missing or defect cable or bad electrodes

More details in the info line are presented in the statistical analysis of treatments in chapter "Control of Therapy" (page 60).

Real Time Clock Setup

Before first use of **PeroBravo** the “Real Time Clock” that is integrated should be set. It’s necessary for control of future therapy sessions. **PeroBravo** has a small battery that will support the clock with energy for several years. For this reason, there is no need to reset the clock, even if the 4 main batteries have been removed.

For setting time and date, please select any custom program 17-20 and press the **SET**-button for at least 5 seconds. You’ll reach the configuration-mode. Now press **SET** once again for 10 seconds. HOUR is indicated and you can set the hour with the right +button (00-23).

Press the left +button to confirm this choice and step forward to setting MINUTE. This is indicated on the screen as well, giving the opportunity to select the minutes via right +button in the same way as before.

This way you are able to set one after the other:

- HOUR
- MINUTE
- DAY
- MONTH and
- YEAR

Press the **PRG**-button on the frontside of the device to save the settings. All values will be automatically stored and the clock starts working.

Language Selection

PeroBravo can be set to one of 6 different languages. If in the setting-mode of the RTC (see above) SET is pressed for 5 seconds once again - the language selection screen is displayed. The following options are available:

- German
- English
- French
- Spanish
- Italian
- Polish

Finally you store your choice with the **PRG**-button.

Preset Programs - Single Pulses

PRG 1-6 contain long single pulses, where the polarity changes after each pulse. This change helps to prevent skin irritation. Here it is normal that every second pulse feels slightly weaker, despite the same intensity. Due to the length of the pulses, please use rubber electrodes in well moistened sponge bags only.

This kind of treatment is particularly intended for initial therapy of paresis to prevent a muscular atrophy in the regeneration phase of the affected nerve. If the activation of the muscle improves during treatment, the success of therapy will not be restricted later by weakness of the muscles.

PRG 1-3:

These programs work with triangular single pulses of different pulse widths. The shape of triangle can be changed each time. It always starts with an ascent of 90% but it is possible to adjust the ascent in steps of 5% between 50% (isosceles) and 90%. The ascent is stored for future uses.

PRG 4-6:

These are trapezoidal single pulses of different lengths. The shape of the trapezoid is variable in each case. An ascent of 80% and a top level of 10% is preset, but it is possible to adjust the ascent from 0% to 80% in increments of 5%. This change only affects the duration of the top-level since the descent of the trapezoids is always 10% of the pulse duration. With increasing treatment success, the ascent can be shortened more and more and thus approximated gradually to a pulse more rectangular in shape. If the ascent is shortened during a running program, the intensity will decrease proportionally for safety reasons and must be adjusted once again if necessary. Each program defaults to 80% ascent each time the program is started.

Any adjustment of the shape of pulses in the programs 1-6 is displayed both visually and in numerical values. All changes of intensity are indicated as well. The increase of intensity between two single pulses is limited to max. 3 mA (see page 38 - Safety Adjustments).

Program No.	Form of the pulses	Pulse width	Pause between the pulses	Number of repetitions	maximum intensity
PRG 01	Triangle	600 ms	5 secs.	120	45 mA
PRG 02	Triangle	300 ms	4 secs.	150	55 mA
PRG 03	Triangle	100 ms	3 secs.	200	60 mA
PRG 04	Trapez	600 ms	5 secs.	120	45 mA
PRG 05	Trapez	300 ms	4 secs.	150	55 mA
PRG 06	Trapez	100 ms	3 secs.	200	60 mA

If the SET button is pressed for 5 seconds or more after selecting the respective program, you can choose automatic or manual triggering of the pulse with the right +button. However, the minimum time before a pulse can be triggered is restricted to the programs pausetime and cannot be changed.

Preset Programs - Biphasic Current

Programs 7-10 are biphasic rectangular pulses with different pulse widths and frequencies.

PRG 7:

This program contains so called LIB-current (according to Mokrusch). LIB means „long impulse bidirectional“. These are quite long-lasting pulses of 40 ms length (positive and negative component) that are delivered with a frequency of 10 Hz for 10 seconds. After a pause of 70 seconds the second of six treatment-phases starts. Due to the large muscular effort this program should not be used more often than twice a day.

PRG 8:

PRG 8 uses biphasic current, which has a pulse length of 1 ms (positive and negative component) and a frequency of 50 Hz. The duration of these pulses is still longer than it is normally used for muscle stimulation of healthy muscles.

PRG 9:

The pulse width used here is even shorter than in PRG 08, and the values are only slightly above the pulse width that is normally used for stimulation of non-denervated muscles.

Program-No.	PRG 07	PRG 08	PRG 09
	LIB-current		
Form of pulses	Rectangle biphasic	Rectangle biphasic	Rectangle biphasic
Pulse width (positive part)	20 ms	500 μ S	400 μ S
Frequency	10 Hz	50 Hz	50 Hz
Work Time	10 secs.	5 secs.	8 secs.
Rest Time	70 secs.	10 secs.	16 secs.
Ramp up	1,2 secs.	0,8 secs.	1 secs.
Ramp down	1,2 secs.	0,8 secs.	1 secs.
Total time of program	8 min.	10 min.	20 min.

Preset Program - Biphasic Current

PRG 10:

Program 10 is a multi-phase stimulation program for healthy muscles - promoting blood circulation for warm-up, then work-/rest-rhythms for muscle training, followed by a purification phase for accelerated regeneration.

Program-No.	PRG 10				
	Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	3 min.	8 min.	3 min.	8 min.	8 min.
Form of pulses	Rectangle biphasic				
Pulse width (positive part)	150 μ S	200 μ S	150 μ S	200 μ S	200 μ S
Frequency	5 Hz	45 Hz	5 Hz	45 Hz	8-12 Hz
Mode	continuous	Work/Rest	continuous	Work/Rest	modulated
Work Time		4 secs.		4 secs.	
Rest Time		10 secs.		10 secs.	
Ramp up		1 secs.		1 secs.	
Ramp down		1 secs.		1 secs.	
Total time of program	30 min.				

Preset Program - Monophasic Current

PRG 11:

Previous studies have shown that good results in the stimulation of peroneal palsy can also be achieved by monophasic pulses. For this reason program 11 offers monophasic positive pulses (according to Vodovnik). Since skin irritation by production of electrolyte can not be excluded, self-adhesive electrodes may not be used, only rubber-electrodes with wet sponge bags. In addition, a very short total duration of the program has been selected and the intensity is limited to 40 mA. If the treatment time is not sufficient for the intended training effect, the program can be restarted.

Program-No.	PRG 11
Form of pulses	Rectangle monophasic positive
Pulse width	20 ms
Frequency	25 Hz
Work Time	3 secs.
Rest Time	4 secs.
Ramp up	0,5 secs.
Ramp down	0,5 secs.
Total time of program	3 min.

Preset Program - Diagnostic

PRG 12:

Program 12 is used to create a strength-duration-curve (SDC). It consists of 24 sequences with monophasic positive triangular and rectangular pulses in 12 different pulse widths each. In the first sequence rectangular pulses are sent at an interval of 2 seconds, the intensity can be increased by pressing the left +button until the first response of the muscle is visible or palpable. The intensity increases in steps of 1 mA. Between two pulses it can be increased by a maximum of 3 mA. The value is added to the coordinate sheet (page 85). Then the second sequence starts again at 1 mA, by pressing the right + / - buttons (%). This second sequence delivers triangular pulses of the same length and the above-described procedure is repeated. In this way, on the sheet of coordinates a triangle-curve is produced step by step as well as a rectangle curve. These curves can then be analysed as specified in the relevant literature. We recommend to make some copies of this page before using it.

Program No.	PRG 12				
	Diagnostic program for creating an I/t-curve				
Form of pulses	Rectangle monophasic positive (RIC) / Triangle monophasic positive (TIC)				
Sequence No.	Pulse Width	Form	Pause	max. Intensity	Display
01	1000 ms	Rectangle	2 secs.	30 mA	RIC
02	1000 ms	Triangle	2 secs.	30 mA	TIC
03	800 ms	Rectangle	2 secs.	35 mA	RIC
04	800 ms	Triangle	2 secs.	35 mA	TIC
05	600 ms	Rectangle	2 secs.	45 mA	RIC
06	600 ms	Triangle	2 secs.	45 mA	TIC
07	500 ms	Rectangle	2 secs.	50 mA	RIC
08	500 ms	Triangle	2 secs.	50 mA	TIC
09	300 ms	Rectangle	2 secs.	60 mA	RIC
10	300 ms	Triangle	2 secs.	60 mA	TIC
11	200 ms	Rectangle	2 secs.	60 mA	RIC
12	200 ms	Triangle	2 secs.	60 mA	TIC
13	100 ms	Rectangle	2 secs.	60 mA	RIC
14	100 ms	Triangle	2 secs.	60 mA	TIC
15	50 ms	Rectangle	2 secs.	60 mA	RIC
16	50 ms	Triangle	2 secs.	60 mA	TIC
17	10 ms	Rectangle	2 secs.	60 mA	RIC
18	10 ms	Triangle	2 secs.	60 mA	TIC
19	5 ms	Rectangle	2 secs.	60 mA	RIC
20	5 ms	Triangle	2 secs.	60 mA	TIC
21	3 ms	Rectangle	2 secs.	60 mA	RIC
22	3 ms	Triangle	2 secs.	60 mA	TIC
23	1 ms	Rectangle	2 secs.	60 mA	RIC
24	1 ms	Triangle	2 secs.	60 mA	TIC

The Programs - Training of Walking

PRG 13-16 are designed for training of walking by current with biphasic rectangular pulses only. All the other programs become deactivated if the foot-switch is connected (s.p. 58).

Besides three pre-set programs (13-15) with different treatment-parameters PRG 16 offers the possibility to set custom parameters according to the therapists requirements.

All programs for training of walking are intended for selective therapeutic practice. The walking exercise is useful for those patients who have already made progress in treating their peroneal palsy or are affected only by a lighter form of foot dorsiflexion. It requires a high degree of concentration on the therapy.

At the beginning of training the foot switch is usually placed in the shoe under the heel of the paretic foot. Only in exceptional cases, when spasticity does not allow the heel strike, the switch should be placed under the footpad of the healthy foot. In this case it is necessary to change the operating-mode from HEEL to TOE. This is done by pressing the SET button for at least 5 seconds after selecting the desired program - then HEEL or TOE can be chosen with the right +button (%).

While in „HEEL“ position the switch triggers current only when released, it is reversed in „TOE“ position, the current flows if the switch is pressed.

In order to achieve a better gait, an adjustment of a „start delay time“ is possible by using the left +button. In this way not only the beginning of the stimulation is delayed but also the end of stimulation. If this is unwanted „stop delay“ can be adjusted separately - without influence on „start delay“

All individual selection or adjustments can be finally confirmed by pressing PRG.

PRG 13-15:

The intensity is adjusted once with the left + /-button to an effective but comfortable level, prior to the start of training. The program starts 5 seconds after the last button press wherein the control of current flow is taken over by the foot-switch. **Intensity cannot be changed during the running program anymore. If an increase or decrease of intensity is wanted, the program must be re-started.**

Program-No.	PRG 13	PRG 14	PRG 15
Form of pulses	Rectangle biphasic	Rectangle biphasic	Rectangle biphasic
Pulse width (positive part)	200 μ S	350 μ S	180 μ S
Frequency	30 Hz	36 Hz	45 Hz
Stim Start Delay / Stim Stop Delay	400 ms (0 ms - 1000 ms adjustable)		
Ramp up	0,4 secs.	0,4 secs.	0,4 secs.
Ramp down	0,4 secs.	0,4 secs.	0,4 secs.
Total time of program	20 min.	20 min.	20 min.

Custom Program - Walking Training

PRG 16:

This program is available only if the foot-switch is connected. The parameters that determine success of treatment can be selected according to individual therapist requirements. All statements about the programs 13-15 shall apply here.

Programming can be started by pressing and holding the SET-button for five seconds. Each parameter that can be altered is indicated by small flashing bars in the display. Now changes can be made with the right +/-buttons (%). If the desired value is adjusted the choice is confirmed with the left +button and the next item that can be changed is indicated.

The following parameters can be selected in this program.

- Frequency of stimulation (Hz)
- Pulse width (μ S)
- Ramp up time of current
- Ramp down time of current
- Start delay
- Stop delay
- Heel- or Toe-mode

Program-No.	PRG 16	Steps
Form of pulses	Rectangle biphasic	
Frequency	10-60 Hz	1 Hz
Pulse width (positive part)	50-500 μ S	10 μ S
Ramp up	0-500 ms	10 ms
Ramp down	0-500 ms	10 ms
Stim Start Delay / Stim Stop Delay	0-1000 ms	10 ms
Mode	Heel /Toe	

Custom Programs - Single Pulses

Programs 17-18 offer the chance to create custom treatment-programs containing single pulses according to therapist requirements. The process of programming is identical to program 16: Programming is started by holding the SET-button for five seconds. Each function that may be altered will be indicated by small flashing bars. Then these values can be changed with the right +/--buttons (%). If the value is adjusted, the choice is confirmed with the left +button and the next parameter is indicated.

PRG 17: Triangular Pulses

PRG 18: Trapezoidal Pulses

The following parameters can be adapted to individual needs here:

- Shape of pulses
- Pulse width (ms)
 - 5 ms - 300 ms = Minimum 3 secs.
 - 301 ms - 600 ms = Minimum 4 secs.
 - 601 ms - 1000 ms = Minimum 5 secs.
- Pause between pulses
- Number of repetitions
- Automatic or manual trigger of pulses

If „manual trigger“ is selected, each pulse is triggered by pressing the SET-button or mA+ or mA-. If mA+/- is used the intensity of the following pulse will increase/decrease by 1 mA.

A pulse cannot be triggered during the pause period.

Press the PRG-button to save this configuration. The program is stored now until it is reconfigured.

PRG 17 + PRG 18		
Parameter	Values	Steps
Pulse width	5 - 1000 ms	1 - 10 ms
Pause between the pulses	3 - 60 secs.	1 sec.
Number of repetitions	5 - 200	1

The use of rubber electrodes with wet sponge bags is required for these programs.

Custom Programs

Biphasic Current - Monophasic Current

PRG 19-20:

These are multiphase custom programs with biphasic rectangular (19) or monophasic positive rectangular pulses (20). When using program 20 please note that its intensity is limited to 40 mA and that the use of rubber electrodes with wet sponge bags is required.

Programming starts when the SET-button is pressed for at least 5 seconds. Small flashing bars indicate the function that can be altered in the display. Now changes can be made by using the right +/-buttons (%). If the desired value is set, the choice is confirmed with the left +button. Now the next variable function is indicated. When a phase is completed, the programming of the next phase can be started by pressing SET. In this way it is possible to program up to 5 phases in different modes. The following functions are available:

- Duration of the phase (min.)
- Mode (Continuous, Burst, Modulated or Work/Rest)
- Frequency (Hz)
- Pulse width (μ S/ms)

If a modulated Mode is active additional settings are required:

- Frequency low limit (Hz)
- Frequency high limit (Hz)
- Pulse width low limit (μ S)
- Pulse width high limit (μ S)
- Intensity modulation (50-100%) - 100% = no modulation of intensity
- Modulation time (secs.)

If a work/rest rhythm is chosen it's necessary to determine the following additional values:

- Work time (secs.)
- Rest time (secs.)
- Ramp up (secs.)
- Ramp down (secs.)

In PRG 19 the programming of the Work/Rest-mode is limited to 40 minutes in total. Each single Work/Rest-phase is limited to max. 20 minutes. These setting restrictions are valid for all modes in PRG 20 as well, while PRG 19 allows unlimited programming of continuous, modulated or Burst-mode. If in any of these two programs the selected pulse width for Work/Rest-phases is setup above 400 μ S, the Work/Rest settings will be limited depending on the safety limits (s.p. 57)

Pressing SET at any time saves the programming of the current phase and selects the next phase. Pressing PRG saves the completed program and returns to home screen.

If any phase time is set to zero minutes, the current phase and all following phases are disabled.

Custom Programs

Biphasic Current - Monophasic Current

PRG 19 + PRG 20		
Parameter	Values	Steps
Phase time	1 - 20 minutes in work/rest mode unlimited in all other modes	1 min.
Mode	Continuous, Burst, Modulated, Work/Rest	
Continuous		
Frequency	1 - 120 Hz	1 - 5 Hz
Pulse width (positive part)	50 - 350 μ s	10 μ s
Burst		
Frequency	35 - 120 Hz	1 Hz
Pulse width (positive part)	50 - 200 μ s	10 μ s
Modulated		
Frequency low	1 - 120 Hz	1 - 5 Hz
Frequency high	1 - 120 Hz	1 - 5 Hz
Pulse width low	50 - 350 μ s	10 μ s
Pulse width high	50 - 350 μ s	10 μ s
Intensity modulation	50 - 100 %	10%
Modulationtime	2 - 60 secs.	1 - 10 secs.
Work/Rest		
Frequency	1 - 120 Hz	1 - 5 Hz
Pulse width (positive part)	50 μ s - 35 ms	10 μ s / 1 ms
Worktime	2 - 30 secs.	1 sec.
Resttime	2 - 90 secs.	1 sec.
Ramp up	0,1 - 10 secs.	0,1 sec.
Ramp down	0,1 - 10 secs.	0,1 sec.

Depending on the selected pulse width (PW) the following safety limits have been established for the protection of the user:

- PW 50 μ S - 400 μ S:
no restrictions
- PW 410 μ S - 5 ms:
maximum frequency = 50 Hz
maximum work time = 10 secs.
minimum rest time = work time x 2 (if PW 410 μ S - 1000 μ S)
minimum rest time = work time x 5 (if PW 2 ms - 35ms)
maximum phase time = 10 min.
- PW 6 ms - 35 ms:
maximum frequency = 250/number of ms
maximum work time = 10 secs.
minimum rest time = work time x 5
maximum phase time = 10 min.

Foot-Switch


The switch that is delivered with the **PeroBravo** is a „Force Sensing Resistor“ (FSR) which realizes whether it is pressed or not. Depending on the setting for heel- (default) or toe-operation the current is turned off or triggered. The connection of the foot-switch to the jack that is named FOOT in the head of the device makes the use of the programs 13-16 possible. All the other programs become deactivated if the foot-switch is connected.

If the foot switch is removed while PRG 13-16 is selected, CONNECT will flash in the info line of the display.

Usually the switch is placed in the shoe under the heel of the paretic foot (see drawing below). Only in the very rare cases that the heel doesn't allow any contact, the switch may be placed under the footpad of the healthy foot. In this case TOE-operation must be adjusted (s.p. 53).

Grey = paretic foot
White = healthy foot

So that the switch does not slip during the training of walking, we recommend securing it with a suitable adhesive tape to the desired position.

The switch is intended for use on one single patient only.



Electrode Placement Paresis



For walking exercises in programs 13-16 it is recommended using the included self-adhesive electrodes, since these are better protected against accidental slipping during movement. These programs contain only short biphasic pulses, in which skin-irritation is minimal.

When using programs with longer pulse widths ($> 500 \mu\text{s}$) we recommend to work with rubber electrodes and wet sponge bags only. This is definitely necessary in the programs 11 and 20 since long monophasic pulses are used.

Regarding placement of electrodes for the treatment of peroneal palsy or for gait-training, two different methods have been proven.

PeroBravo was also designed for the treatment of the following forms of paresis.

The device can also be used for the treatment of facial paralysis. Since the bony anchorage of the muscles in the face is only weak there is an associated risk of over exertion. For this reason we recommend to use only low current intensities and to place the electrodes only in accordance with the detailed instructions of an experienced physician. A general placement recommendation is not given here.

Control of Therapy

The **PeroBravo** records all treatment sessions and collects statistics. In this way doctors or therapists have the chance to get to know the possible reasons for the presence or the absence of treatment success. Longer sections of unintended use can be detected and prevented. The device stores a total of 365 days per 5 sessions. When 365 days of data are reached, consecutive days will overwrite the earliest days. The day number will continue increasing (366, 367, ...). For this reason the statistics will always show the actual data of the last 365 treatment days.

To read the memory, the UNLOCK-button in the battery-compartment (s.p. 45) must be pressed twice - if the device has been locked to a single program - or once, if the device has been unlocked. Now **STATS. GL** (for „global statistics“) will be displayed in the info line. If the device has not been used at all, the screen will show **NO STATS**.

Global Statistics

Global statistics show the average values of all recorded sessions. The actual term is named in the info line and two small bars are flashing beside the respective value to make assignment easier. The change from one term to the next is done with the right +button (%). In this way the following information can be read one after the other:

- the most used PROGRAM
- the average used intensity in mA (AVG. MA)
- the number of treatment days (DAYS.USED)
- the number of treatment hours in total (HRS. USED)
- the average number of sessions per day (AVG. SESS)

By pressing the left +button (mA) the more detailed daily statistics are available. The screen shows **STATS.DAY** and the information regarding the first session of the first day can be seen.

Pressing left +button (mA) = next day / Pressing left -button (mA) = previous day

Pressing SET = next session of the selected day

Pressing right +/-buttons = scrolling through the session details.

Session statistics

- number of treatment day (STATS.DAY)
- average intensity used in this session (AVG. MA)
- length of this session (SESS. LEN)
- used program (PROGRAM)
- date of this session (SESS.DATE)
- start time of this session (SESS.TIME)

Holding the left + or - button (mA) pressed for three seconds will return to „Global Statistics“. If the SET-button is pressed for 5 seconds, **DEL.STATS** will be displayed and after confirmation with SET - for 5 seconds once again - all statistics will be deleted and all individual programmed parameters are set to their default values.

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To achieve the best possible effect, one of the electrodes should be put on the motor point of the respective muscle if possible. Try to find out this point by slightly moving the electrode around. The other electrode is fastened on the remote muscle end. However, you also obtain a good effect if you put the electrodes respectively on the opposite muscle ends. Please make sure that you always use rubber electrodes with wet sponge bags with monophasic pulses.

■ = red cable connection and ■ = black cable connection

While stimulating, please keep this device in reach, to be able to finish the treatment at once, if current becomes uncomfortable.

Please follow the references of this manual.

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

In most cases it is very easy to find the suitable points for the placing of the electrodes. Placing the electrodes either side of the painful area gives good results. We nevertheless recommend to coordinate the points of electrode-placement with the doctor or physiotherapists if possible. The distance between the electrodes shouldn't be much more then 20 cm. The polarity of the electrodes does not have any significance in TENS therapy. Please make sure that you always use rubber electrodes with wet sponge bags with monophasic pulses. You see some examples of possible placement-points on the following pages.

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














Batteries and Battery-Charger

PeroBravo is supplied with rechargeable batteries. If these batteries don't have sufficient charge to operate, you can recharge using the charger. For this purpose, please insert four rechargeable NiMH-batteries of type AA with right direction for polarity, as it is drawn in the charger and then plug it to the mains (230V). Charging starts immediately and after about 5 hours, the batteries are fully charged. There is no danger of overcharging because the device has a corresponding protection feature. While charging is in progress, the small yellow LED is blinking slowly. If this LED is continuously on, the batteries are charged completely. Now disconnect the unit from the power supply and then remove the charged batteries and put them into the battery compartment of **PeroBravo**.

If the yellow LED is flashing quickly, the charger stops charging because a wrong or damaged battery has been detected and needs to be removed.

If you want to use other than the supplied batteries, please make sure that they have a minimum capacity of 1800 mAh or more - to avoid frequent charging.

Safety Instructions:

- Don't charge batteries of different size (AA or AAA) at the same time.
- Please use NiMH-batteries only and do not try to charge normal batteries (primary batteries) - this could cause an explosion.
- Always charge 2 or 4 batteries, never 1 or 3 batteries. If you only load 2 batteries, please set them next to each other without a gap in the left or right charging slots, and not in the middle.
- Avoid the simultaneous charging of old and new batteries or charging of batteries with different capacity or of different manufacturers.
- Please never open batteries by force and don't throw them into a fire.
- The device may only be operated in enclosed, dry areas - it has to be protected from rain and moisture.
- It is normal that the batteries become hot when charging.
- If the power plug or the cabinet show any damage, the charger may not be operated.
- The charger may not be opened or altered, repairs shall be done by qualified persons only.
- Do not operate the charger unsupervised and disconnect it from the mains after use. Keep it out of reach of children or persons with reduced physical or mental capacity.
- Cleaning and maintenance may only be done after disconnection.

The address or website of the manufacturer can be found on the rear label of the charger. The actual charging current depends on the type and number of batteries to be charged. For loading the supplied four AA batteries it is 525 mA.

Please see also the instruction manual of the manufacturer.

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In accordance with European regulations [1], used electrical and electronic equipment may not be discarded in household or unsorted municipal waste. Use separate collection facilities. Contact your local government or council for information regarding the collection systems available. If electrical appliances are

disposed of in landfills or dumps, hazardous substances can leak into the ground water and get into the food chain damaging the health and well being of people and animals.

[1] European Community directive 2002/96/EC on waste electrical and electronic equipment.

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The batteries should be recharged or replaced when they are drained or stop functioning. The batteries supplied with the **PeroBravo** may be recharged many times. If these batteries have reached the end of their lifetime or if non-rechargeable batteries are used please remove them from the unit when they are empty and follow your local recycling guidelines.

Literatur / Literature

Barrett, C.L.; Mann, G.L.; Taylor, P.N. and Strike, P. – A randomized trial to investigate the effects of functional electrical stimulation and therapeutic exercise on walking performance for people with multiple sclerosis; *Multiple Sclerosis* 2009; 15: 493–504

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Breen, Paul P.; Corley, Gavin J.; O’Keeffe, Derek T.; Conway, Richard; ÓLaighin, Gearóid - A programmable and portable NMES device for drop foot correction and blood flow assist applications, *Medical Engineering & Physics* 31 (2009) 400–408

Breen, P.P.; O’Keeffe, D.T.; Conway, R.; Lyons, G.M. - A system for the delivery of programmable, adaptive stimulation intensity envelopes for drop foot correction applications; *Medical Engineering & Physics* 28 (2006) 177–186

Byrne, C.A.; O’Keeffe, D.T.; Donnelly, A.E.; Lyons; G.M. - Effect of walking speed changes on tibialis anterior EMG during healthy gait for FES envelope design in drop foot correction; *Journal of Electromyography and Kinesiology* 17 (2007) 605–616

Crochetiere, W.J.; Vodovnik, L.; Reswick, J.B. - Electrical stimulation of skeletal muscle - a study of muscle as an actuator; *Med. & biol. Engng.* Vol. 5, pp. 111-125, 1967.

Hausdorff, J.M.; Ring, H. - Effects of a New Radio Frequency– Controlled Neuroprosthesis on Gait Symmetry and Rhythmicity in Patients With Chronic Hemiparesis; *Am J Phys Med Rehabil* 2008; 87: 4–13

Hui-Chan, C.W.Y.; Ng, S.S.M.; Mak, M.K.Y. – Effectiveness of a home-based rehabilitation programme on lower limb functions after stroke; *Hong Kong Med J* Vol 15 No 3 Suppl. 4, 2009

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Lyons, Gerard M.; Sinkjær, Thomas; Burridge, Jane H.; Wilcox, David J. - A Review of Portable FES-Based Neural Orthoses for the Correction of Drop Foot; *IEEE Transactions on Neural Systems and Rehabilitation Engineering*, Vol. 10, No. 4, December 2002

Maležič, M.; Bogataj, U.; Gros, N.; Dečman, I.; Vrtačnik, P.; Kljajić, M.; Aćimović-Janežič, R. - Application of a programmable dual-channel adaptive electrical stimulation system for the control and analysis of gait; *Journal of Rehabilitation Research and Development*, Vol. 29 No. 4, 1992, 41 -53

Mansfield, Avril; Lyons, Gerard M. - The use of accelerometry to detect heel contact events for use as a sensor in FES assisted walking; *Medical Engineering & Physics* 25 (2003) 879–885

Mokrusch, T. - LIB-Stimulation (Long Impulse, Bidirectional) Improves Reinnervation in Patients With Long Term Denervation; In: Karenik, T.; Veltink, P.; Jaeger, R.; editors;

Literatur / Literature

Proceedings of the 7th Annual Conference of the International Functional Electrical Stimulation Society; Ljubljana, Slovenia, 2002: 144-6

Ng, S.S.M.; Hui-Chan, C.W.Y. - Does the use of TENS increase the effectiveness of exercise for improving walking after stroke? A randomized controlled clinical trial; *Clinical Rehabilitation* 2009; 23: 1093–1103

O’Keeffe, Derek T.; Donnelly, Alan E.; Lyons, Gerard M. - The Development of a Potential Optimized Stimulation Intensity Envelope for Drop Foot Applications; *IEEE Transactions on Neural Systems and Rehabilitation Engineering*, Vol. 11, No. 3, September 2003

Pappas I.P.I.; Keller, T.; Mangold, S.; Popovic, M.R.; Dietz, V.; Morari, M. - A Reliable Gyroscope-Based Gait-Phase Detection Sensor Embedded in a Shoe Insole; *IEEE Sensors Journal*, Vol. 4, No. 2, April 2004

Sabut, S. K.; Kumar, R.; Mahadevappa, M. - Design of a programmable multi-pattern FES system for restoring foot drop in stroke rehabilitation, *Journal of Medical Engineering & Technology*, Vol. 34, No. 3, April 2010, 217–223

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Sheffler, L.R.; Hennessey, M.T.; Naples, G.G.; Chae, J. - Peroneal Nerve Stimulation versus an Ankle Foot Orthosis for Correction of Footdrop in Stroke: Impact on Functional Ambulation; *Neurorehabilitation and Neural Repair* 20(3); 2006

Singh, B.; Xu, Q.-G.; Franz, C.K.; Zhang, R.; Dalton, C.; Gordon, T.; Verge, V.M.K.; Midha, R.; Zochodne, D.W. - Accelerated axon outgrowth, guidance and target reinnervation across nerve transection gaps following a brief electrical stimulation paradigm; *J Neurosurg*, Vol. 116; 2012; 498-512

Stein, R.B.; Everaert, D.G.; Thompson, A.K.; Chong, S.L.; Whittaker, M.; Robertson, J.; Kuether, G. - Long-Term Therapeutic and Orthotic Effects of a Foot Drop Stimulator on Walking Performance in Progressive and Nonprogressive Neurological Disorders; *Neurorehabilitation and Neural Repair*; 24(2); 152–167

van Swigchem, R.; Vloothuis, J.; den Boer, J.; Weerdesteyn, V.; Geurts, A.C.H. – Is transcutaneous peroneal stimulation beneficial to patients with chronic stroke using an ankle-foot orthosis? A within-subjects study of patients’ satisfaction, walking speed and physical activity level; *J Rehabil Med*, 2010; 42: 117–121

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Electromagnetic Compatibility

Guidance and manufacturer's declaration - electromagnetic emissions		
The PeroBravo unit is intended for use in the electromagnetic environment specified below. The customer or the user of the PeroBravo should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The PeroBravo uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The PeroBravo is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not Applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable	

Guidance and manufacturer's declaration – electromagnetic immunity			
The PeroBravo unit is intended for use in the electromagnetic environment specified below. The customer or the user of the PeroBravo should ensure that it is used in such an environment.			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/-6 kV contact +/-8 kV air	+/-6 kV contact +/-8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Electromagnetic Compatibility

Guidance and manufacturer's declaration – electromagnetic immunity

The **PeroBravo** unit is intended for use in the electromagnetic environment specified below. The customer or the user of the **PeroBravo** should ensure that it is used in such an environment.

Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the PeroBravo unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ 80MHz to 800MHz $d=2.3 \sqrt{P}$ 800MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: <div style="text-align: right; margin-top: 10px;">  </div>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m 80 MHz to 2,5 GHz	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

A) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **PeroBravo** unit is used exceeds the applicable RF compliance level above, the **PeroBravo** unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the **PeroBravo** unit.

B) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Electromagnetic Compatibility

Recommended separation distances between portable and mobile RF communications equipment and the PeroBravo unit.

The **PeroBravo** unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the **PeroBravo** unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **PeroBravo** unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[\frac{3,5}{V_1}\right]\sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3,5}{E_1}\right]\sqrt{P}$	800 MHz to 2,5 GHz $d = \left[\frac{7}{E_1}\right]\sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Reserved for changes.


Printing or duplication of this manual - also in an abridged version - is forbidden and can only be allowed with a written approval of “MTR+ Vertriebs GmbH”.

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