Neuromuskulärer Stimulator
Neuromuscular Stimulator

RehaBravo

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Content

Please read the Operators Manual carefully before using the device
and pay attention to all references !

	Page
I. Neuromuscular Stimulator RehaBravo	
Delivery & Accessories	33
Introduction - What is RehaBravo ?	34
Indication - Contraindication	35
Warnings	36
Advice on Application of Electrodes	37
Technical Data - Explanation of Symbols	38
Maintenance - Warranty	39
Setup Real Time Clock - Safety Adjustments	40
II. Display Elements of the RehaBravo	
Display and Operating Elements	41
Description of the Device	42
Battery Compartment with LOCK-button	43
III. Application of the RehaBravo	
Preset Programs	44-45
Setting the Custom Programs	46-47
IV. Muscle Stimulation	
Types of Muscle Fibres	48-49
Effect of Training to the Muscular Fibres	50
Frequency Selection	51
Programs 1 - 15	52-56
Abbreviations	57
V. Control of Therapy	
Lock- & Unlock-Buttons - Statistics	58
Reading the Statistics	59-60
Placement of Electrodes	
Electrode-Placement for Muscle-Stimulation	61-70
Electrode-Placement for TENS-Therapy	71-75
Literature	76-77
Waste Disposal - Recycling of Batteries	78
Anatomy	79
Electromagnetic Compatibility	83-85
Copyright	85

Delivery & Accessories

Included in the delivery:

REF 010E-002: 1 device **RehaBravo**, 2 lead wires, 4 self adhesive electrodes (40 x 80 mm), 4 batteries type AA and 1 operators manual

The expected service life of **RehaBravo** is 10 years. Service life of accessories is one year.

Accessories:

REF

011E-006: 1 set self adhesive electrodes 40x40 (4 pieces)

011E-007: 1 set self adhesive electrodes 40x80 (4 pieces)

011E-051: 1 set self adhesive electrodes 32mm Ø (4 pieces)

011E-008: 1 vaginal probe “IncoProbe V”

011E-009: 1 rectal probe “IncoProbe R”

011E-012: 2 lead wires

011F-016: 4 batteries type AA

011E-046: 1 operators manual **RehaBravo**

011E-050: 1 carrying case

The batteries life is approximately 12 hours or more at 20 mA per channel. The storage time of the batteries is approximately 2 years. The service life (expiry date) of electrodes is printed on each package.

RehaBravo shall only be used with original accessories.

Introduction – What is RehaBravo ?

RehaBravo is a neuromuscular stimulator designed and developed by MTR+ in close working relationship with medics and patients, taking most recent knowledge and developments into consideration. This device is a modern 2 channel stimulator, which allows simultaneous use of 2 pairs of electrodes (= 4 electrodes) at one time. If necessary, you can use vaginal- or rectal-probes as well. **RehaBravo** provides a highly effective muscle stimulation for therapeutic purposes.

The special emphasis of **RehaBravo** is the orientation to the field of orthopaedic rehabilitation. Based on years of experience in the medical field, special value was placed on careful and gentle muscle build up as well as effective regeneration by relaxation, purge and active recovery besides the numerous programs which the device provides for the increase of muscle power and fatigue resistance.

It also contains a special program for pain therapy. The “Dense-Disperse-Stimulation”, created by Prof. Han offers a permanent alternation of frequency and pulse width in a 3-second rhythm.

Additionally the three freely configurable programs enable the specialist to compile specially tailored programs for their own individual purpose. This characteristic of the **RehaBravo** offers each physician or physiotherapist, and also to the experienced patient, the possibility of a focused coordinated treatment for the respective state of health. As these programs, can be changed at any time, for example to be altered according to the progress of therapy, an almost unrestricted variation range is available with this device to the user.

The small size of this neat but powerful unit makes transport easier. **RehaBravo** fits into any bag and can be used in holiday as well as at home without any assistance.

After termination of the therapy, **RehaBravo** can be re-used by other persons. Please clean it before with a mild disinfectant cleaner.

The clip on the back of the unit allows to fasten the device at the clothes but the unit is intended to be supported by hand during normal use. The belt clip is for temporary relief only.

Indications

RehaBravo was designed for the following indications:

Muscle Stimulation:

- Avoidance and delay of muscular spasms (PRG 1, 2, 13)
- Relaxation of muscular spasms (PRG 2, 13)
- Avoidance and delay of muscular amyotrophy (PRG 9, 10)
- Recovery of the muscle (PRG 3, 4, 5, 6, 9, 10, 11, 12)
- Conserve and extend the mobility (PRG 3, 4, 5, 6)
- Extension of the capillary bed (PRG 1, 2, 14)
- Speeding up the process of removal of toxic substances from metabolism (PRG 1, 2, 13, 14)
- Strengthening of the muscle (PRG 3, 4, 5, 6, 9, 10, 11, 12)
- Activation of former inactive muscular fibres (PRG 7, 8)
- Increase of the power of resistance of the muscle against tiredness (PRG 3)
- Avoidance of muscular imbalances (PRG 3, 4, 5, 6, 9, 10, 11, 12)

Nerve Stimulation:

- Increase of release of bodily endorphins
- Relieve or elimination of pain, for instance in the case of:
 - Cervical Spine pain
 - Thoracic spine pain
 - Lumbar spine pain
 - Sporting accidents
 - Circulation problems
 - Neuralgia
 - Myalgia
 - Scar and phantom pain
 - Arthritis
 - Fracture pain

For all these indications you can use program 15.

Contraindications

RehaBravo may not be used in the following cases:

- By patients with electronic implants like pumps or a pace maker because functional disturbances can occur
- By patients with myocard damages (damage to the cardiac muscle) or arrhythmia (inconsistent heart beat)
- In presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.
- When operating machinery with a potential of danger
- By patients with epilepsy
- When inflammation of the vagina, anus or urinary tract is present
- By patients with mechanical ileus
- By patients with febrile illnesses or infectious diseases

Warnings

Please pay unconditional attention to the following basic rules:

- Do not stimulate transcerebral (via the head).
- Do not stimulate via the Carotis-Sinus-Nerve (carotid artery).
- Do not stimulate in the area of the Sinus knot (heart).
- Do not stimulate on the Bulbous (eye).
- Do not stimulate in the Larynx-Pharynx-Area (larynx /pharynx).
- Do not stimulate in places with skin defects.
- Do not position the electrodes directly on the spine.

Take special care with the use of the electric stimulation:

- if after an operation the healing process can be interrupted through a muscle contraction.
- after acute bleeding or bone fractures with a tendency to bleed.
- during menstruation or pregnancy
- in case of desensitized skin with nerve damage.
- if stimulation is located close to osteosynthesis material (metal in the body).
- if patients cannot use the device properly, due to mental or physical limitations.

Please note necessarily the following information:

RehaBravo is a type BF equipment and may not be immersed in liquids.

If at the same time the connection is made with a surgical high frequency apparatus, this could result in burn marks underneath the electrodes.

If device operation occurs close to short waves or a microwave, instability of the initial values of the stimulator could occur.

The application of the electrodes nearby the thorax can increase the risk of cardiac fibrillation. Should skin irritation occur whilst using the **RehaBravo**, please interrupt the application and inform your doctor - to exclude an allergic reaction.

If in the case of treatment with the probe, pain or irritation appear in the area of the pelvis, or bleeding in-between menstruation, interrupt the therapy and consult your doctor.

Patients with metal implants may only stimulate with permission of an experienced doctor.

Stimulation in the area of tumors is allowed only by instruction of the specialized physician.

The **RehaBravo** may not be operated if the battery compartment is not closed.

Before opening the battery compartment, make sure that there are no cables in the jacks at the head of the unit to eliminate any possibility of current flow.

The operator should not touch the battery contacts & patient simultaneously.

Never connect the device directly to a battery charger or to any other mains powered equipment.

If using rechargeable batteries be sure to use a CE approved battery charger.

The **RehaBravo** is protected against debris larger than 12.5 mm in diameter from entering the cabinet, is not water resistant and rated IP20.

If you have swallowed any part of the device or any of its accessories, please contact a doctor immediately.

Electrodes can contain amounts of nickel. Allergy sufferers should check this before using.

Please keep this device stored out of the reach of children.

Advice on Application of Electrodes

Make sure that **RehaBravo** is turned off (no indication visible on the LCD display).

Now connect the single poled ends of the lead wires with the electrodes.

The output plugs of the electrode cables will be plugged into the channels A or B of the **RehaBravo**.

Furthermore make sure, that all plugs fit completely into the electrode sockets as well as the device sockets.

If you are using a probe, lubricate it with water or a gel suitable for electrotherapy.

Please don't use any oil, ointment or cream, since all these have an isolating effect.

Now place the electrodes or the probe according to the instructions of your doctor or physiotherapist.

In the case of the skin gel losing its adhesion, you can't achieve more than 6 mA - please replace the electrodes. The same effect, however, also occurs when using a damaged cable.

Therefore, if the normal treatment is not possible by using new electrodes, the cable may be faulty.

If you use a probe, please clean it after use carefully with a normal detergent and dry it completely.

Electrodes/probes which were damaged due to an improper handling may not be used any more.

Electrodes/probes with a damaged cable isolation (e.g. visible copper wire) may not be used any more.

Electrodes with damaged, missing or stained aluminium coating or carbon layer, also with insufficient adhesive gel may not be used.

Electrodes/probes may not be altered (e.g. by means of scissors or other tools).

Do not pull the electrodes/probes by force on the electrode cable.

Apply the electrodes on clean and healthy skin only.

In the case of strong bodily hair growth, the hair should be shortened with scissors.

Neither stretch the skin nor the electrodes when applying the electrodes.

To remove the electrodes from the skin, please lift the electrodes at the rim towards the middle.

Effective current densities of more than 2 mA/cm² require an increased attention of the user with regard to the intensity adjustment.

If you use self-adhesive electrodes, we recommend the use of our self-adhesive electrodes with an area of 40x80 mm.

In case of pelvic floor treatment we recommend the use of our vaginal- or rectal-probes. If you should use a probe, please clean it carefully before and after use and dry it completely.

Technical Data

Channels: Two isolated independent channels

Waveform: Asymmetrical, rectangular, biphasic with zero DC-current

Intensity: 0 - 90 mA into 500 ohm load,

70 mA maximum into 1000 ohm load

65 mA maximum into 1500 ohm load

(Indication only, actual intensity varies depending on condition of the electrodes)

Classification: Internally powered, constant current, continuous operation, type BF applied part, class IIa (MDD 93/42/EEC), IP 20

Frequency: 1 - 200 Hz

Pulse Width: 50 μ s - 450 μ s

Ramp up/down: 0,1 - 9,9 seconds

Time adjustment: 1 minute - 9 hours, 59 minutes

Power supply: 4 batteries 1,5 V type AA

Voltage max.: 180 V without load (80 V into 1000 ohm load)

Dimensions: L 141,6 mm x W 74 mm x H 40,4 mm

Weight: 250 g (incl. batteries)

Environmental conditions for storage & transportation: 0-93% atmospheric humidity,

-25 to +70 degree centigrade. Environmental conditions for use: 0-93% atmospheric humidity,

+5 to +40 degree centigrade, air pressure 700-1060 hPa.

This device is quality controlled and fulfills the rules and regulations set by the German law for medical products (MPG) as well as the guidelines of the European Community 93/42/EEC for medical products. Therefore this device carries the CE-sign „CE 0123“. The serial number is in the battery compartment.

Manufacturer: MTR+ Vertriebs GmbH,
Kamenzer Damm 78, 12249 Berlin - www.mtrplus.com

Explanation of Symbols

Attention
Warning

Type BF applied parts

Serial-Number
of the product

Keep dry

Batch code

Manufacturer

Date of
manufacture

CE-mark: conforms to
essential requirements of
the MDD 93/42/EEC

Not suited for persons
with cardiac pacemaker

Read the instructions

Maintenance

RehaBravo has been designed to be maintenance free. Nevertheless, please adhere to the following advices:

If needed, please clean **RehaBravo** with a mild detergent, water and a damp cloth. Do not expose **RehaBravo** to extreme humidity and do not submerge under water. You may clean the cables with a damp cloth. You may also from time to time treat the cables with talcum powder in order to prevent cracks and to further on going usage. A clean cloth, moistened with an alcohol free disinfection liquid can be used to clean the device. Please follow all instructions of the manufacturer of the liquid.

After termination of the therapy, **RehaBravo** can be re-used by other patients. Please clean the device with a mild disinfectant cleaner.

The **RehaBravo** does not require calibration. The **RehaBravo** is not repairable and contains no user serveable parts. No modification to the **RehaBravo** is allowed.

When changing the batteries, please make sure to insert the new batteries within 30 minutes after removal of the spent batteries. We recommend that you always have spare batteries with the device.

Warranty

MTR+ Vertriebs GmbH provides a 24 month warranty, commencing from the date of our delivery (invoice date), that in case of failure, the **RehaBravo** will be serviced or replaced free of charge. All returns must first be authorised by MTR+ Vertriebs GmbH in advance.

Excluded from this warranty are failures of the **RehaBravo** which result from excessive use, careless or improper handling as well as damage by force.

The warranty loses its lawful meaning if the producer MTR+ Vertriebs GmbH has not given its consent in writing prior to repairworks or servicing being carried out. Without this written consent the manufacturer will not return the device.

The manufacturer reserves the right of technical changes, which may serve to increase treatment success, the functionality of the device or its accessories.

Setup Real Time Clock

RehaBravo has a “Real Time Clock”, that shows year, month, day, hours and minutes to allow the therapeutic memory the temporal assignment of detailed treatment data.

For setting time and date, please select any custom program 16-18 and press the **SET**-button for at least 5 seconds. You'll reach the configuration-mode. Now press **SET** once again for 10 seconds.

With the +/- button of channel A (left side) you can step forward to the following values

- hour (indicated by: H)
- minute (indicated by: MI)
- day (indicated by: dA)
- month (indicated by: MO)
- year (indicated by: Y)

The respective value is flashing and can be changed with the +/- buttons of channel B (right side). Finally you save all your settings by pressing the PRG-button once again.

Before starting the therapy, please set the clock once if necessary.

If **RehaBravo** has been without batteries for more than 30 minutes, the clock must be set again.

Safety Adjustments

1. The +buttons become inactive after one minute running time of a program to avoid an unintentional increase of intensity. This lock is switched off again, only by pressing the +button three-times repeatedly. Now an increasing of intensity is possible again for another minute.
2. If the electrical circuit should be interrupted during operation (> 6 mA), for example by losing an electrode or by a broken lead wire, the device stops stimulating immediately.
3. If any phase comes to its end, the intensity falls down to 8 mA automatically and must be reset by the user to the desired setting.

Display and Operating Elements

„ON/OFF“: By pressing the button once, the device will either be switched on or off.

„+/- buttons“: This button adjusts the intensity for the channels A and B. The device also has two independently controllable output channels A and B. When pressing the respective arrow buttons the intensity will either be increased or decreased in steps of 1 mA. 60 seconds after the last intensity regulation, a further increase may not be carried out due to an electronic barrier. Pressing the „+“ button three times, will release the lock and the intensity can be increased again. A reduction of the intensity is possible at any time.

Whilst configuring the programs 16 – 18 one switches between the programmable parameters by pressing the left „+/-“ button (channel A), and through pressing the right „+/-“ button (channel B) the desired value/data for respective parameter can be set. In this way, the application time, the stimulation, the pause, rising pulse and falling pulse, as well as synchronous or asynchronous stimulation can be selected.

„PRG“: With this button the preset programs 01–15 or the three custom programs 16–18 can be selected for individual use.

Pressing the PRG-button once at any time during a running program will reduce the mA on both channels to zero and pause the clock display with the actual treatment time. The clock continues to show the time of further treatment but it is frozen with the “:” symbol no longer flashing. In addition, this button stores the current adjustment in the programming mode and ends it. Then the device returns to the home screen of the selected program.


„SET“: With the SET-button you can select the different program phases. With each press of this button, the phase in progress is terminated and the next program phase is started. This applies to either the running of the programs 01-15 and also to the programming and the sequence of programs 16–18. Programs with one single phase will be finished by pressing the SET-button.

If you keep the SET-button pressed for 5 seconds whilst being in programs 16–18, the symbol of the clock is flashing and you will be able to start programming of your individual settings. Pressing SET once again for 10 seconds allows to setup the RTC.


Holding the „SET“ button for 5 seconds while the statistics are displayed, will display „dEL“. Holding „SET“ a second time for 5 seconds will delete the statistics, set all custom programs to default values and show the RTC setup screen.

Description of the Device


Channel A **Channel B**
Places for Lead Wires

Battery Compartment with LOCK-button

RehaBravo allows the doctor or therapist to lock one program for the patient to use (s. p. 58). In this way the medical practitioner can be quite sure that the patient only uses the device with the correct program. The LOCK-button can be pressed with a pin of the cable or any other suited thin object. Pressing will fix the selected program so that different programs cannot be used anymore.

UNLOCK button removes the LOCK of programs and shows global statistics

Arrangement of the batteries

The serial number is located on the back of the unit above the rear label.

Preset Programs 01–15

Preparations:

Make sure that **RehaBravo** is switched off. Cleanse the area of the skin which is to be treated with clean water. Following that, make sure that all plugs are completely connected to the respective electrode sockets as well as device sockets. Place the the adhesive electrodes or the vaginal-/rectal-probe according to the advices in pages 35-37 of this manual.


RehaBravo is operated by 4 batteries of type AA. Open the battery compartment at the back by pushing down the lid and insert the batteries. Please ensure the correct polarity in accordance with the markings. **RehaBravo** is now operational. If the battery voltage is too low, the battery-symbol in the LCD-display will flash. Please replace the battery at this time. Empty batteries constitute special waste and therefore have to be disposed of at an appropriate collection point. Never dispose of batteries via domestic or household waste.

Adjustment:

1. Press the button „ON/OFF” once.

RehaBravo is now switched on. The last used program appears on the LCD-display.


Parameters programmed for the first phase will be displayed as well as the quantity of all phases, the actual treatment mode and the total treatment time.

Preset Programs 01–15

2. By pressing the program button „PRG” you can choose from the programs 01 – 15 (pre-programmed).
3. This is followed by pressing the buttons „+/-” below the channel used.
Now the selected program runs automatically and the LCD-screen changes.

The current phase will be indicated by an flashing arrow and in the upper right corner you can see the remaining phase time. In the bottom line the intensity used in each channel appears. If the button „+” is pressed repeatedly the intensity of stimulation increases. If the button „-” is pressed repeatedly the intensity of stimulation is reduced. For the setting of the intensity it is vital that one feels the electric current but that discomfort is never incurred.

Please take into account that for safety reasons the intensity falls back to 8 mA at the beginning of each new phase and must be increased again.

Each phase end is indicated by 5 peeptones.

4. If you wish to interrupt the program, then press the button „PRG” once. The program will run again for the remaining program time when pressing the button „+”. In this case, the intensity must be reset to the desired level.
5. By pressing the „SET” button you end the current phase and move to the next. Programs with one single phase will be finished.
6. **RehaBravo** is switched off when pressing the button „On/Off” or automatically four minutes after the program is finished.

Setting the Custom Programs 16-18

1. Switch on **RehaBravo** with „ON/OFF” button.
2. Choose the desired program (16-18) with „PRG” button.
3. Press „SET” button for at least five seconds. **RehaBravo** changes to the programming mode and the upper watch, which indicates phase duration, flashes to show that this function can be altered.

In the following please use the left „+” button always then, if you want to store your choice and to change to the next step. The right „+/-” buttons are respectively used in order to alter the chosen and flashing parameter.

4. Select with the right „+/-”buttons the desired time for the first phase.
5. By pressing the left „+”button you will be able to store your choice and change to the next step. In the display the mode now flashes, which you adjust according to your wishes by pushing the right „+/-”buttons.

If you have chosen „CONT” or „BURST” you can determine frequency and after storage and change with the left +button respectively the pulse width by pushing the right „+/-” buttons.

If you have chosen „WORK/REST”- rhythm you can determine the following parameters one after the other by pushing the right „+/-”buttons.

- frequency (**Hz**)
- pulse width (**µS**)
- work time (**sec**)
- rest time (**sec** = recovery time)
- increase / decrease ramp up time of the pulse (**RU**)
- increase / decrease ramp down time of the pulse (**Rd**)
- synchronous (**SY**) or alternating (**AL**) output of the channels A and B. If you have selected “**SY**” you can program **delay-time (dL)** of 0 sec. up to 5 sec. from channel A to B.

If you have chosen „MOD” you can alter the following parameters one after the other.

- low limit of frequency (**Hz LO**)
- upper limit of frequency (**Hz HI**)
- low limit of pulse width (**µS LO**)
- upper limit of pulse width (**µS HI**)
- low limit of intensity (AdJ) - the display in the right corner below shows values from 0.5 to 1.0. This adjusts the intensity of current within the modulation time between 50% (=0.5) and 100% (=1.0) of the original set mA-value. The intensity, of a modulated program, can never increase above the 100% of the set value.
- modulation time (possible choice: 2 sec. up to 60 sec)

Setting the Custom Programs 16-18

6. Push the „**SET**” button to select the second phase of the program. You can program up to 5 phases of your own in each of the three configurable programs. If you want to reduce a configured program of 5 phases to less phases, you simply set the time of the first unwanted phase to zero.
7. The programming is completed by pressing the „**PRG**” button and the program is stored.

Selectable frequency and pulse width settings:

Continuous-Mode: 1 to 200 Hz, 50 to 450 μ S

Burst-Mode: 35 to 200 Hz, 50 to 200 μ S

Modulated-Mode: 1 to 200 Hz, 50 to 450 μ S

Work/Rest-Mode: 1 to 200 Hz, 50 to 450 μ S

It is not possible to program more than 300 μ S if the frequency is above 100 Hz.

The settings can be made as follows.

- Frequency: 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 15 Hz, and then increasing in 5 Hz steps until max. 200 Hz
- Pulse width: 50 to 450 μ S in steps of 10 microseconds
- Intensity: 50-100% in steps of 10%
- Modulation time: 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 15, 20, 30, 40, 50, 60 seconds
- Work time: 2-99 seconds in steps of 1 second
- Ramp up time: 0,1-9,9 seconds in steps of 1 second
- Ramp down time: 0,1-9,9 seconds in steps of 1 second
- Rest time: 2-99 seconds in steps of 1 second
- Delay: 0,0-5,0 seconds in steps of 0,1 second
- Session time: 1-599 minutes in steps of 1 minute

Muscular Stimulation

It is a known fact that nerves control muscles by transmission of neurological codes. This code or message, depending on the required type of muscular fibre, will be transmitted in varying frequency bands. The muscle begins to contract when it receives electrical stimulus, whether the impulse is generated by the brain or by electrical stimulation. The physiological method of neuromuscular stimulation requires impulses, that are similar in condition to the naturally occurring nerve signals. Provided that the imitation of nature is as exact as possible, electrical stimulation may, if necessary, be applied over a long periods of time, without causing side effects. Muscular system of the skeleton (horizontal striped musculus) consist of numerous long and thin fibres, muscular fibres that span between tendons, by which they are linked to the bones. The respective suitable stimulus is dependant on the type of muscular fibre, that is to be reached. The main distinctions between the different fibre type is as follows:

Type I (red muscular system)

This type of fibre is also known as ST-fibre (slow twitch fibres – slowly contracting fibres). The necessary gain of energy for the functioning of the muscles happens aerobically, that is through oxidative metabolism. The motoneuron that stimulates this fibre, has a slow transmitting speed. These type of fibres are thin and of red colour (which refers back to the existence of myoglobin molecular). Internally these fibres contain a high number of mitochondria and oxidative enzymes. Type I muscular system is extremely resistant to tiredness, as it is responsible for all types of activities of tonic nature, slow and are connected with the support of bodily posture. These slow fibres are surrounded by a dense capillary network, which allows optimal functioning of the aerobic metabolism by prolonged activity in connection with a negligible expenditure of energy. The red muscle fibres bestow stability on the bodily posture muscles and support the joint. This muscle system is very important for all endurance sports like cycling, long distance running, swimming etc.

Type IIa (white muscle system)

This type is also known as FOG-fibres (fast twitch oxidative glycolytic fibres – fast contracting fibres with oxidative glycolytic metabolism). These are stimulated by a motoneuron of the phased type that has a higher transmitting speed than tonic motoneuron. The fibres are white because of the lack of myoglobin and have a mixed metabolic activity. They are rich in glycogen and glycolytic enzymes and also contain mitochondria enzymes, the whole metabolism is more anaerobic than the aerobe-oxidative one. These fibres are also equipped with a capillary network, that transports the necessary oxygen for the aerobic process. Fibres of type IIa are therefore capable of fast contractions, which are marked by a significant expenditure of energy, that should also be maintained for long periods of time. The fibres are therefore relatively resistant to tiredness.

Muscular Stimulation

Type IIb (white muscular system)

This type is also known as FG-fibres (fast twitch glycol fibres – fast contracting fibres with glycolytic metabolism). This fibre is stimulated by a phasal motoneuron with a cellular body and a very big axon that transmits impulses into the muscle at very high velocity. These fibres are white and contain a high level of glycogen and glycolytic enzymes, that enables them to release a high level of anaerobic energy. The contraction is quite fast and produces great power. Due to lack of mitochondrion, these fibres will not be able to endure long-lasting activities and fatigue easily, especially in an untrained muscle. Fibres of type IIb play a very important role in all human activities, which require an explosive expenditure of energy, e.g. sports like sprinting, weight lifting, swimming, jumping etc.

Spread of muscular fibres

Nearly all human muscles show a mixed form of different fibre types. The aforementioned types of fibres appear in diverse percentages in the muscles. The relation between fibres of type I and type II can vary considerably. The muscles for posture (back and stomach muscular systems) show a higher proportion of type I fibres, whereas muscles for movement / motion (leg and arm muscular systems) usually consist mainly of type II fibres.

Effect of Training to the Muscular Fibres

Enhancement of physical capabilities by training or electrical stimulation can be related back to the fact, that the human body is able to adjust to a certain extent to changing environmental conditions. If the stress to a muscle is increased beyond a threshold value, a corresponding adaptation process will ensue. The increase of muscular power is initially due to improved intro-muscular co-ordination, which means that for one movement more fibres will be activated simultaneously than before. Only in the second step will the muscular mass be increased. Whether this occurs exclusively due to hypertrophy of the muscular fibres (meaning the enlargement of the fibre diameter) or if it is the result of hyperplasie (meaning the increase of muscular fibres) is scientifically not proven. It is probable that both phenomena concurrently are responsible for the muscular growth.

The aforementioned training effects to the muscular system can be attained in a more passive way by electrical stimulation. One should not make the mistake of neglecting the physical training in favor of electrical stimulation. A movement, as the aim of almost every sporting activity, consists of a finely tuned response of different muscles or groups of muscles. Smooth movements will only be achieved through an optimum of coordination of contracting the active muscles (agonists) while relaxing their opponents (antagonists). These smooth movements are unconditional requirements for sporting success. The necessary neuromuscular coordination cannot be trained with this kind of muscle stimulation.

To avoid muscular imbalances, it is urgently advisable to apply strengthening programs for both the agonists and the antagonists in equivalent amounts.

Beside the strengthening of muscles, the support of muscular blood circulation is an additional main target of electrical stimulation. A higher capillary density enables the improved blood circulation. Hence an increase supply of nutritive substances, especially oxygen, is made available and the aerobic metabolism benefits. Electrical stimulation not only allows an improvement of the supply, but also supports the fast disposal of metabolic waste products like lactic acid or carbon dioxide. Speeding up physical regeneration is of vital importance regarding the efficiency of athletes - especially in phases of hard training or competition.

A very short burst of electrical stimulation produces only a short contraction or „single shock“ after which the muscle falls back into its natural form and length, which it possesses in the resting position. If the stimulation occurs sequentially, we realise that the contracting effects caused by overlapping of the contracting phases are additive, as the mechanical twitches last distinctly longer than the electrical stimulation. This phenomena is called „incomplete tetanus“. Neither „single shock“ nor „incomplete tetanus“ will normally be observed in a voluntary human movement.

The phenomena of „complete tetanus“ means the state of a muscular contraction which is caused by repeated electrical stimulation of the motoric nerves with a sufficient frequency. Hereby the single pulses merge and are no longer distinguishable. In this situation, the muscle contracts - a measurable power is exercised at the tendon endings.

Frequency Selection

5 Hz or lower Single twitches can be achieved with this frequency, after which the muscle can relax completely for a short time. These frequencies are used in the starting phase of the strengthening programs with the purpose of preparing the muscle for the coming working phase without tiring it. A frequency of 3 Hz lies within the frequency band for the production of bodily endorphins for the relief of pain and general relaxation.

5 – 15 Hz These frequencies are also referred to as shaking frequencies, which will not lead to the „complete tetanus“ but otherwise do not allow for complete relaxation in-between the impulses. This frequency range is selected for improvement of muscular tonus, support and stability of the joints.

10 Hz is the natural frequency of the slow oxidative muscle fibre (Type I). The electrical stimulation increases the resistance to tiredness by improving the density of the capillary system and it increases the ability of muscles to tolerate an oxygen-reduction. Onerous metabolic waste will be further reduced through the „milking of the veins“. This frequency band is of great importance in the regeneration phase after high levels of stress and can be used daily over a number of hours for sportive or therapeutic purpose.

15 – 30 Hz This frequency can be used for increase of stamina of the muscle. This is the natural frequency band of the oxidative-glycolytic muscle fibres (Type IIa). The frequency can be applied up to one hour daily.

30 – 60 Hz The muscle stays in a state of complete tetanus from a frequency of 30 Hz upwards. These frequencies are applied for the strengthening of the muscle and recruitment of the fast glycolytic muscle fibres (Type IIb).

A treatment in these frequency bands should only be applied for relatively short periods because the tiring of the muscle will begin within a few minutes of electrical stimulation.

60 – 120 Hz These frequencies will often be selected in order to train the muscle to produce great energy within a very short time (explosive force). This kind of stimulation should only be applied for very short periods.

Programs

Progr.: 01		Warming up				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	6				
Mode		Cont				
Frequency work	Hz	5				
Frequency rest	Hz					
Pulse duration	μ S	300				
Modulation time	secs					
Ramp up time	secs					
Ramp down time	secs					
Work time	secs					
Rest time	secs					
Overall time	6 min					

Progr.: 02		Capillary				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	20				
Mode		Cont				
Frequency work	Hz	10				
Frequency rest	Hz					
Pulse duration	μ S	250				
Modulation time	secs					
Ramp up time	secs					
Ramp down time	secs					
Work time	secs					
Rest time	secs					
Overall time	20 min					

Progr.: 03		Increase of muscles 1				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	20				
Mode		W/R				
Frequency work	Hz	50				
Frequency rest	Hz					
Pulse duration	μ S	300				
Modulation time	secs					
Ramp up time	secs	1,2				
Ramp down time	secs	1,2				
Work time	secs	10				
Rest time	secs	20				
Overall time	20 min					

Programs

Progr.: 04		Increase of muscles 2				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	5	15	10		
Mode		Cont	W/R	MOD		
Frequency work	Hz	5	60	8-12		
Frequency rest	Hz		3			
Pulse duration	μ S	300	350	250		
Modulation time	secs			6		
Ramp up time	secs		2			
Ramp down time	secs		1,5			
Work time	secs		7			
Rest time	secs		14			
Overall time	30 min					

Progr.: 05		Strengthening of muscles 1				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	20				
Mode		W/R - MOD				
Frequency work	Hz	40-60				
Frequency rest	Hz	3-5				
Pulse duration	μ S	300				
Modulation time	secs	10				
Ramp up time	secs	1				
Ramp down time	secs	1				
Work time	secs	10				
Rest time	secs	15				
Overall time	20 min					

Progr.: 06		Strengthening of muscles 2				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	5	20	10		
Mode		Cont	W/R	MOD		
Frequency work	Hz	5	75	8-12		
Frequency rest	Hz		3			
Pulse duration	μ S	300	300	250		
Modulation time	secs			6		
Ramp up time	secs		1,5			
Ramp down time	secs		1			
Work time	secs		5			
Rest time	secs		12			
Overall time	35 min					

Programs

Progr.: 07		Severe partial peripheral paralysis				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	10				
Mode		W/R				
Frequency work	Hz	60				
Frequency rest	Hz					
Pulse duration	μ S	450				
Modulation time	secs					
Ramp up time	secs	0,6				
Ramp down time	secs	1				
Work time	secs	4				
Rest time	secs	20				
Overall time	10 min					

Progr.: 08		Medium partial peripheral paralysis				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	15				
Mode		W/R				
Frequency work	Hz	50				
Frequency rest	Hz					
Pulse duration	μ S	300				
Modulation time	secs					
Ramp up time	secs	0,5				
Ramp down time	secs	1				
Work time	secs	5				
Rest time	secs	15				
Overall time	15 min					

Progr.: 09		Severe muscle atrophy				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	10				
Mode		W/R				
Frequency work	Hz	6				
Frequency rest	Hz					
Pulse duration	μ S	350				
Modulation time	secs					
Ramp up time	secs	0,8				
Ramp down time	secs	1				
Work time	secs	5				
Rest time	secs	35				
Overall time	10 min					

Programs

Progr.: 10		Medium muscle atrophy				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	20				
Mode		W/R				
Frequency work	Hz	15				
Frequency rest	Hz					
Pulse duration	µS	300				
Modulation time	secs					
Ramp up time	secs	0,6				
Ramp down time	secs	1				
Work time	secs	6				
Rest time	secs	30				
Overall time	20 min					

Progr.: 11		Scoliosis / Kyphosis				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	60				
Mode		W/R				
Frequency work	Hz	35				
Frequency rest	Hz					
Pulse duration	µS	250				
Modulation time	secs					
Ramp up time	secs	1,5				
Ramp down time	secs	1,5				
Work time	secs	8				
Rest time	secs	20				
Overall time	60 min					

Progr.: 12		Scoliosis 2				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	599				
Mode		WR - MOD				
Frequency work	Hz	20-50				
Frequency rest	Hz					
Pulse duration	µS	250				
Modulation time	secs	60				
Ramp up time	secs	1,5				
Ramp down time	secs	1,5				
Work time	secs	8				
Rest time	secs	20				
Overall time	599 min					

Programs

Progr.: 13		Muscle recovery				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	20	20	10		
Mode		MOD	MOD	W/R		
Frequency work	Hz	2-5	5-10	10		
Frequency rest	Hz					
Pulse duration	µS	250-150	250-150	200		
Modulation time	secs	10	10			
Ramp up time	secs			2		
Ramp down time	secs			2		
Work time	secs			10		
Rest time	secs			10		
Overall time	50 min					

Progr.: 14		Decongestion				
		Phase 1	Phase 2	Phase 3	Phase 4	Phase 5
Phase time	min	20				
Mode		WR - MOD				
Frequency work	Hz	40				
Frequency rest	Hz	2-5				
Pulse duration	µS	250				
Modulation time	secs	10				
Ramp up time	secs	1				
Ramp down time	secs	1				
Work time	secs	5				
Rest time	secs	30				
Overall time	20 min					

Progr.: 15		Pain therapy				
		Phase 1				
Phase time	min	30				
Mode		HAN				
Frequency	Hz	100 Hz / 200 µs changing continuously with				
Pulse duration	µS	2 Hz / 300 µs in 3:3-seconds-rhythm				
Overall time	30 min					

Abbreviations

AdJ:	Adjustment of intensity
AL:	Alternate function of the channels A and B
BURST:	The burst mode delivers 2 blocks of 9 impulses each per second
CONT:	In the continuous mode a permanent unchanged current is used
dA:	Day
dEL	Delete
dL:	Delay between channels A and B
HAN:	Change of frequency and pulse width in a 3:3 seconds rhythm (indicated as CONT + BURST)
HI:	Upper value (high)
Hz:	Hertz is the measurement unit for the frequency of impulses per second.
LO:	Lower value
mA:	Milli-amp is the measurement unit for the electrical current intensity.
MOD:	In the modulated mode the frequency and/or the pulse width and/or the intensity changes in a predefined rhythm automatically.
Rd:	Ramp down
RU:	Ramp up
SEC:	Seconds
SY:	Synchronous function of the channels A and B
TI:	Time
TM:	Time of modulation-interval
W/R:	Work/Rest mode
µS:	Microsecond is the measurement unit for the length of every single impulse.

LOCK-/UNLOCK Buttons

The **RehaBravo** offers the ability to lock one program in place so that other programs or settings cannot be selected or changed. This allows the doctor or therapist to ensure that the patient uses only the program that has been selected for him.

The LOCK switch in the battery compartment (s. p. 43) can be used by means of the pin of the electrode cable. To do this, please disconnect all cables from the device, switch on the device and open the battery compartment. After pressing LOCK you hear a single beep and the selected program is locked. In this case, you cannot use any other program anymore until UNLOCK is pressed.

Pressing the UNLOCK button in the battery compartment (s. p. 43) releases the program lock. A second press shows the global statistics on the display. If the unit was not locked, these statistics are already displayed on the first press.

Statistics

The therapist can read information about the compliance of the patient on the display at any time he wants. This way it is possible to evaluate whether the device was used according to the instructions of the doctor (quantity, intensity, etc.)

RehaBravo offers the chance to get a quick survey about the stored data of all treatment-sessions. The global statistics display the following averages:

- number of treatments per day
- number of days in use
- most used program
- total time of use
- intensity per channel

If more detailed information is required, **RehaBravo** allows the therapist to have a closer look at each single session of the last 60 treatment-days. You can see:


- number and date of the treatment-day
- number of session at this day
- time of session-start
- number of program used in this session
- total-time of the session
- intensity used for each channel in this session


Reading the Statistics

If you want to read the memory, switch on the device and press the UNLOCK-button in the battery compartment twice if the unit was locked in one program or once if it was not locked. The global statistics can be viewed on screen now.

Global Statistics:

To get more detailed information press the + button of the left channel (A). You will see the data of the first session of the first day.

Reading the Statistics

The next sessions of the selected day are displayed - one after the other - by pressing the „SET“-button. Each next day can be obtained by pressing the + button of channel A again. If this button is held for at least 3 seconds, the device returns to global statistics.

If you require information about time and date of one single session please press any of the right +/--buttons (channel B). With these buttons you can scroll between session data, date and time.

PHASE:1 2 3 4 5

Number of chosen session flashes — 03

Treatment day — 09 — Year

— dA for Date

Day — 16 — Month

PHASE:1 2 3 4 5

Number of chosen session flashes — 03

Treatment day — 71 — TI for Time

Hour — 15 — Minute

The data stored in memory is retained after viewing, so that therapy-control is possible anytime. A maximum of 60 days of data can be stored with 5 sessions per day. When 60 days of data are reached, consecutive days will overwrite the earliest days. The day number will continue increasing (61, 62, ...) and will be displayed when the unit is unlocked until the **RehaBravo** is reset.

If you press the PRG button for three seconds, you return to the home screen of the selected program. The device will not be locked. If the therapy shall continue with only one program available, the LOCK-function must be reactivated after reading (s. p. 58).

The memory should be emptied completely before a new person is using the device. After pressing the UNLOCK-button you have to press SET twice for 5 seconds each. All individual settings in programs 16-18 will be reset as well. The **RehaBravo** is set to it's default state and all statistics are deleted.

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To achieve the best possible effect, one of the electrodes should be put on the motor point of the respective muscle if possible. Try to find out this point by slightly moving the electrode around. The other electrode is fastened on the remote muscle end. However, you also obtain a good effect if you put the electrodes respectively on the opposite muscle ends.

■ = red cable connection and ■ = black cable connection

While stimulating, please keep this device in reach, to be able to finish the treatment at once, if current becomes uncomfortable.

Please follow in each case absolutely the references on page 35-37 of this manual.

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In most cases it is very easy to find the suitable points for the placing of the electrodes. It is frequently meaningful to take the hurting place to the middle between the two electrodes. We nevertheless recommend to coordinate the points of electrode-placement with the doctor or physiotherapists if possible. The distance between the electrodes shouldn't be much more then 20 cm. The polarity of the electrodes is without any significance in TENS therapy.

You see some examples of possible placement-points on the following pages.

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









In accordance with European regulations [1], used electrical and electronic equipment may not be discarded in household or unsorted municipal waste. Use separate collection facilities. Contact your local government or council for information regarding the collection systems available. If electrical

appliances are disposed of in landfills or dumps, hazardous substances can leak into the ground water and get into the food chain damaging the health and well being of people and animals.

[1] European Community directive 2002/96/EC on waste electrical and electronic equipment.

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
The batteries should be replaced when they are drained or stop functioning. The batteries supplied with the **RehaBravo** are recyclable. Remove the batteries from the unit and follow your local recycling guidelines.

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Electromagnetic Compatibility

Guidance and manufacturer's declaration - electromagnetic emissions		
<p>The RehaBravo unit is intended for use in the electromagnetic environment specified below. The customer or the user of the RehaBravo should ensure that it is used in such an environment.</p>		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The RehaBravo uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The RehaBravo is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not Applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable	

Guidance and manufacturer's declaration – electromagnetic immunity			
<p>The RehaBravo unit is intended for use in the electromagnetic environment specified below. The customer or the user of the RehaBravo should ensure that it is used in such an environment.</p>			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/-6 kV contact +/-8 kV air	+/-6 kV contact +/-8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Electromagnetic Compatibility

Guidance and manufacturer's declaration – electromagnetic immunity

The **RehaBravo** unit is intended for use in the electromagnetic environment specified below. The customer or the user of the **RehaBravo** should ensure that it is used in such an environment.

Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC 61000-4-6</p>	<p>3 Vrms 150 kHz to 80 MHz</p>	<p>3 Vrms 150 kHz to 80 MHz</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the RehaBravo unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> <p>$d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ 80MHz to 800MHz $d=2.3 \sqrt{P}$ 800MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>
<p>Radiated RF IEC 61000-4-3</p>	<p>3 V/m 80 MHz to 2,5 GHz</p>	<p>3 V/m 80 MHz to 2,5 GHz</p>	<p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: right;">  </div>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

A) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **RehaBravo** unit is used exceeds the applicable RF compliance level above, the **RehaBravo** unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the **RehaBravo** unit.

B) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Electromagnetic Compatibility

Recommended separation distances between portable and mobile RF communications equipment and the RehaBravo unit.

The **RehaBravo** unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the **RehaBravo** unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **RehaBravo** unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
	$d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$	$d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$	$d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Reserved for changes.

Printing or copy of this manual - also in an abridged version - is forbidden and can only be allowed with a written approval of “MTR+ Vertriebs GmbH”.

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