Neuromuskulärer Stimulator  
Neuromuscular Stimulator

# Rehalito

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## Content

Please read the Operators Manual carefully  
before using the device and pay attention  
to all references !

	<b>Page</b>
<b>I. Musclemulator Rehalito</b>	
Delivery & Accessories	22
Indication / Contraindication	23
Warnings	24
Advice on Application of Electrodes	25
Technical Data	26
Explanation of Symbols / Maintenance	27
Warranty / Safety Adjustments	28
Real Time Clock Setup	29
<b>II. Display Elements of the Rehalito</b>	
Display and Operating Elements	30
Description of the Device	31
<b>III. Application of the Rehalito</b>	
Preset Programs	32-33
Setting the Custom Program	34-35
<b>IV. Programs 01-05</b>	36
Abbreviations	37
<b>V. Lock-Button and Statistics</b>	38-39
<b>Placement of Electrodes</b>	40-54
<b>References / Literature</b>	55-56
<b>Waste Disposal / Copyright</b>	57
<b>Anatomy</b>	58
<b>Electromagnetic Compatibility</b>	62-64
<b>Copyright</b>	65

## Delivery & Accessories

### Included in the delivery:

REF: 010E-201: 1 device **Rehalito**, 2 lead wires, 4 self adhesive electrodes (40 x 80 mm), one 9V battery and 1 operators manual

The expected service life of **Rehalito** is 10 years.

### Accessories:

REF

011E-006: 1 set self adhesive electrodes 40x40 (4 pieces)

011E-007: 1 set self adhesive electrodes 40x80 (4 pieces)

011E-051: 1 set self adhesive electrodes 32mm Ø (4 pieces)

011E-008: 1 vaginal probe “IncoProbe V”

011E-009: 1 rectal probe “IncoProbe R”

011F-012: 2 lead wires

011F-013: 2 pieces 9V batteries

011F-037: 1 tube electrode jelly for probes

011E-209: 1 operators manual **Rehalito**

The battery life is approximately 12 hours or more at 20 mA per channel.

The service life (expiry date) of electrodes, probes and electrode jelly is printed on each package.

**Rehalito** shall only be used with original accessories.

## Indications

**Rehalito** was designed for the following indications and is easy to use without any assistance:

### **Muscle Stimulation:**

- Avoidance and delay of muscular spasms (PRG 3)
- Relaxation of muscular spasms (PRG 3)
- Avoidance and delay of muscular amyotrophy (PRG 1, 2)
- Recovery of the muscle (PRG 1, 2)
- Conserve and extend the mobility (PRG 1, 2)
- Extension of the capillary bed (PRG 1, 2, 3, 4)
- Speeding up the process of removal of toxic substances from metabolism (PRG 3)
- Strengthening of the muscle (PRG 1, 2)
- Activation of former inactive muscular fibres (PRG 1, 2)
- Increase resistance of the muscle against fatigue (PRG 1)
- Avoidance of muscular imbalances (PRG 1, 2)

### **Nerve Stimulation:**

- Increase release of bodily endorphins
- Relief or elimination of pain, for instance in the case of:
  - Cervical Spine pain
  - Thoracic spine pain
  - Lumbar spine pain
  - Sporting accidents
  - Circulation problems
  - Neuralgia
  - Myalgia
  - Scar and phantom pain
  - Arthritis
  - Fracture pain

For all these indications you can use program 04.

## Contraindications

### **Rehalito may not be used in the following cases:**

- Patients with a pace maker as functional disturbances can occur
- The **Rehalito** is not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.
- Patients with myocard damages (damage to the cardiac muscle) or arrhythmia (inconsistent heart beat)
- When operating machinery with a potential of danger
- Patients with epilepsy
- When inflammation of the vagina, anus or urinary tract is present (probes)
- Patients with febrile illnesses or infectious diseases

## Warnings

### **Please pay unconditional attention to the following basic rules:**

- Do not stimulate transcerebral (via the head)
- Do not stimulate via the Carotis-Sinus-Nerve (carotid artery)
- Do not stimulate in the area of the Sinus knot (heart)
- Do not stimulate on the Bulbous (eye)
- Do not stimulate in the Larynx-Pharynx-Area (mouth and neck)
- Do not stimulate in places with skin defects or broken skin
- Do not position the electrodes directly on the spine

### **Take special care with the use of the electric stimulation:**

- If after an operation the healing process can be interrupted through a muscle contraction
- After acute bleeding or bone fractures with a tendency to bleed
- During menstruation
- In case of desensitized skin with nerve damage

### **Please note also the following instructions:**

If at the same time the connection is made with a surgical high frequency apparatus, this could result in burn marks underneath the electrodes.

If device operation occurs close to short waves or a microwave, instability of the initial values of the stimulator could occur.

The application of the electrodes nearby the thorax can increase the risk of cardiac fibrillation.

In case of pregnancy, stimulation should only be applied after consultation with a doctor.

If in the case of treatment with the probe, pain or irritation appear in the area of the pelvis, or bleeding in-between menstruation, interrupt the therapy and consult your doctor.

Should skin irritation occur whilst using the **Rehalito**, please interrupt the application and inform your doctor - to exclude an allergic reaction.

Patients with metal and/or electronic implants may only stimulate with permission of an experienced doctor.

The **Rehalito** may not be operated if the battery compartment is not closed.

Before opening the battery compartment, make sure that there are no cables in the jacks at the head of the unit to eliminate any possibility of current flow.

Never connect the device directly to a battery charger or to any other mains powered equipment.

The operator should not touch the battery contacts & patient simultaneously.

The **Rehalito** is protected against debris larger than 12.5 mm in diameter from entering the cabinet and is not water resistant and rated IP20.

**Please, keep this device stored out of the reach of children.**

## Advice on Application of Electrodes

Make sure, that the **Rehalito** is switched off (no display visible on the LCD-Display). Connect the single poled ends of the lead wires to the electrodes plugs with the electrodes. The output sockets of the electrode cables will be plugged into the channels A or B of the **Rehalito**.

Furthermore make sure, that all plugs fit completely into the electrode sockets as well as the device sockets.

Thereafter place the electrodes onto the skin. An adhesive skin gel has been applied to the electrodes, which can be re-used after removal from the skin. In the case of the skin gel loosing its adhesion and you can't achieve more than 6 mA - please replace the electrodes.

The same effect, however, also occurs when using a damaged cable. Therefore, if normal treatment is not possible by using new electrodes, the cable maybe faulty.

Electrodes which are damaged due to careless handling may not be used.

Electrodes with damaged cable insulation (e.g. visible copper wire) may not be used.

Electrodes with damaged, missing or stained aluminium coating or carbon layer, also with insufficient adhesive gel may not be used.

Electrodes may not be altered or tampered with (e.g. by scissors or other tools).

Do not pull the electrodes by force on the electrode cable. To remove the electrodes from the skin, please lift the electrode at the rim towards the middle.

Do not use electrodes with ointments or creams (especially which enhance the blood flow or act as local anaesthesia).

Apply the electrodes on clean and healthy skin only.

In the case of bodily hair growth, the hair should be shortened with scissors.

Neither stretch the skin nor the electrodes when applying the electrodes.

Effective current densities of more than 2 mA/cm<sup>2</sup> require an increased attention of the user with regard to the intensity adjustment.

For standard treatment we recommend the use of our self-adhesive electrodes with an area of 40x80 mm.

If you should use a probe, please clean it carefully before and after use and dry it completely.

## Technical Data

Channels:	Two isolated independent channels
Waveform:	Asymmetrical, rectangular, bi-phasic with zero DC current
Intensity:	0 - 80 mA into 500 ohm load 70 mA maximum into 1000 Ohm load 65 mA maximum into 1500 Ohm load (Indication only, actual intensity varies depending on condition of skin electrodes)
Classification:	Internally powered, constant current, continuous operation, Type BF applied part, class IIa (MDD 93/42/EEC), IP20
Frequency:	1 - 150 Hz (Burst 35 - 150 Hz)
Pulse Width:	50 - 450 $\mu$ S (Burst 50 - 200 $\mu$ S)
Time Settings:	Variable 1 minute to 9 hours 59 minutes (custom program)
Max. Voltage:	180 V without load (70V with 1000 ohm load)
Battery;	9 Volt 6LR61 / 6AM6 recommended

Environmental conditions for storage & transportation:

0-93% atmospheric humidity, -25 to +70 degree centigrade.

Environmental conditions for use: 0-93% atmospheric humidity, +5 to +40 degree centigrade, air pressure 700-1060 hPa.

Dimensions of the device: H 11,7 cm / W 6,6 cm / D 3,3 cm

Weight of the unit: 150 g (incl. battery)

This device is quality controlled and fulfils the rules and regulations set by the German law for medical products (MPG) as well as the guidelines of the European Community 93/42/EEC for medical products. Therefore this device carries the CE-sign „CE 0123“. The serial number is on the rear of the device.

**Manufacturer:** MTR+ Vertriebs GmbH,  
Kamenzer Damm 78, 12249 Berlin - [www.mtrplus.com](http://www.mtrplus.com)

## Explanation of Symbols

Attention  
Warning

Type BF applied parts

Serial-Number  
of the product

Keep dry

LOT number

Manufacturer

Date of  
manufacture

CE-mark: conforms to  
essential requirements of  
the MDD 93/42/EEC

Not suited for persons  
with cardiac pacemaker

Read the instructions

## Maintenance

**Rehalito has been designed to be maintenance free. Nevertheless, please adhere to the following advice:**

If needed, please clean **Rehalito** with a mild detergent, water and a damp cloth. Do not expose **Rehalito** to extreme humidity and do not submerge under water. You may clean the cables with a damp cloth. You may also from time to time treat the cables with talcum powder in order to prevent cracks and to further ongoing usage.

A clean cloth, moistened with an alcohol free disinfection liquid can be used to clean the device. Please follow all instructions of the manufacturer of the liquid.

After termination of the therapy, **Rehalito** can be re-used by other persons. Please clean it before with a mild disinfectant cleaner.

When changing the battery, please make sure to insert the new batteries within 30 minutes after removal of the empty battery. We recommend that you always have a spare battery with the device.

The **Rehalito** does not require calibration.

The **Rehalito** is not repairable and contains no user serviceable parts.

No modification to the **Rehalito** is allowed.

## Warranty

MTR+ Vertriebs GmbH provides a 24 month warranty, commencing from the date of our delivery (invoice date), that in case of failure, the **Rehalito** will be serviced or replaced free of charge. All returns must first be authorised by MTR+ Vertriebs GmbH in advance.


Excluded from this warranty are failures of the **Rehalito** which result from excessive use, careless or improper handling as well as damage by force.

The warranty loses its lawful meaning if the producer MTR+ Vertriebs GmbH has not given its consent in writing prior to repairworks or servicing being carried out. Without this written consent the manufacturer will not return the device.

The manufacturer reserves the right of technical changes, which may serve to increase treatment success, the functionality of the device or its accessories.

## Safety Adjustments

1. The +buttons become inactive after one minute running time of a program to avoid an unintentional increase of intensity . This lock is switched off again by pressing the +button three-times repeatedly. An increasing of intensity is possible again for another minute.
2. If the electrical circuit should be interrupted during operation ( $> 6$  mA), for example by losing an electrode or by a broken lead wire, the device stops stimulating immediately and the display shows the flashing symbol of an electrode. The effected channel will also flash.
3. If any phase comes to its end the intensity falls to 8 mA automatically and must be increased by the user to the desired intensity once again.

## Real Time Clock Setup

**Rehalito** has a “Real Time Clock”, that shows year, month, day, hours and minutes and allows to store detailed statistics of treatment data.

To set time and date, please select the custom program 05 and press the **S**-button for at least 5 seconds. You'll reach the configuration-mode. Now press **S** once again for 10 seconds.

With the +/- button of channel A (left side) you can step forward to the following values

- hour (indicated by: H)
- minute (indicated by: MI)
- day (indicated by: dA)
- month (indicated by: MO)
- year (indicated by: Y)

The respective value is flashing and can be changed with the +/- buttons of channel B (right side). Finally you save all your settings by pressing the PRG-button once again.

Before starting the therapy, please set the clock once if necessary.

If **Rehalito** has been without battery for more than 30 minutes, the clock may need to be reset.

## Display and Operating Elements

**„On/Off“:** By pressing the button once, the device will either be switched on or off.

**„+/- buttons“:** These buttons adjust the intensity for channels A and B. The device has two independently controllable output channels A and B. When pressing the buttons the intensity will either be increased or decreased in steps of 1 mA. 60 seconds after the last change of intensity, a further increase may not be carried out due to an electronic barrier. Pressing the „+“ button three times releases the barrier and the intensity can be altered again.

A reduction of the intensity is possible at any time.

Whilst configuring program PC05 one switches between the programmable parameters by pressing the left „+/-“ button (channel A), and through pressing the right „+/-“ button (channel B) the desired value/data for respective parameter can be set. This way, the application time, the stimulation-mode, the frequency and the pulse width can be selected.

**„P“:** This button selects preset programs 01–04 or custom program PC05.

Press the „P“-button once at any time during a running program to reduce the mA on both channels to zero and pause the program. Pressing „P“ a second time ends the program

Holding „P“ for three seconds while statistics are displayed will exit the statistics screen without deleting them. **Rehalito** will return to the home screen of the selected program and the lock will be cancelled.

The „P“ button stores all adjustments made to the custom program mode and returns to the home screen.

**„S“:** This button is used to setup custom program 05. Press the „S“-button for at least 3 seconds in program PC05 to enter the custom program setting. Pressing „S“ a second time for 10 seconds will enter the real time clock (RTC) setup.

Press the „S“ button once at any time during a running program, to reduce the mA on both channels to zero and pause the program. Pressing „S“ a second time will resume the program and the current will increase over a period of one second back to the previously set intensity level.



Holding the „S“ button for 5 seconds while the statistics are displayed, will display „dEL“.

Holding „S“ a second time for 5 seconds will delete the statistics, set the custom program to default values and show the RTC setup screen. Turn off the unit if you do not want to delete the statistics.

The unit is intended to be supported by hand during normal use. The belt clip is for temporary relief only.

## Description of the Device

## Preset Programs (01 - 04)

### Preparations:

Make sure that **Rehalito** is switched off. Cleanse the area of the skin which is to be treated with clean water. Following that, make sure that all plugs are completely connected to the respective electrode sockets as well as device sockets. Place the vaginal-/rectal-probe or the adhesive electrodes onto the area that has been cleansed, according to the advices in pages 23-25 of this manual.


**Rehalito** is operated by a 9 volt block battery. Open the battery compartment at the back by pushing down the lid and insert the battery by hand. Please ensure the correct polarity in accordance with the markings. **Rehalito** is now operational.

If the battery voltage is too low, the battery-symbol in the LCD-display will flash. Please replace the battery at this time. Empty batteries are special waste and therefore have to be disposed off at an appropriate collection point or recycled. Never dispose of batteries via domestic or household waste.

### Adjustment:

1. Press the button „On/Off” once.


**Rehalito** is now switched on. The program used last, appears on the LCD-display. Program parameters will be displayed as well as the treatment mode and the total treatment time.

2. By pressing the program button „P” you can choose from the programs 01 – 04 (pre-programmed).
3. This is followed by pressing the buttons „+/-” below the channel used. The selected program runs automatically and the LCD displays the intensity as milliamps (mA).

In the upper row you can see the remaining program time and in the bottom line the intensity being applied to each channel. If the button „+” is pressed repeatedly the intensity of stimulation increases. If the button „-” is pressed repeatedly the intensity of stimulation is reduced. For the setting of the intensity it is vital that one feels the electric current but that discomfort never occurs.

4. If you wish to interrupt the program, press the button „S” once and the program will pause. The program will resume for the remaining program time when pressing this button once again. In this case, the intensity increases within one second to the previous used value.
5. By pressing the „P”-button you pause as well, you can continue with new intensity adjustment by pressing the +button. If you press „P” a second time the program is finished.
6. **Rehalito** is switched off when pressing the button „On/Off” or automatically one minute after the program is finished.

## Setting the Custom Program (05)

1. Switch on **Rehalito** by pressing the „**On/Off**” button.
2. Choose the program PC05 with the „**P**” button.
3. Press the „**S**” button for at least three seconds. **Rehalito** changes to the programming mode and the time that indicates program duration, flashes to show that this function can be altered.

**In the following please use the left „+“ button to store your choice and to change to the next step. The right „+/-“ buttons are respectively used in order to alter the chosen and flashing parameter.**

4. Select with the right „+/-”buttons the desired time for the phase.
5. By pressing the left „+”button you will be able to store your choice and change to the next step. Mode now flashes on the LCD and you are able to select the mode that you require by pressing the right „+/-”buttons.

If you have chosen „**CONT**” or „**BURST**” you can set the frequency (**Hz**) and pulse width (**μS**).

If you have chosen „**MOD**” you can set the following parameters:

- upper limit of frequency (**Hz HI**)
- low limit of frequency (**Hz LO**)
- upper limit of pulse width (**μS HI**)
- low limit of pulse width (**μS LO**)
- low limit of intensity (**Adj**) - the display in the right corner below shows values from 0.5 to 1.0. This means that the intensity of current can increase and decrease within the modulation time between 50% (=0.5) and 100% (=1.0) of the originally set mA-value. The upper limit of the intensity (100%) is always the mA value adjusted by the user originally. The intensity in a modulated program, can never increase above the 100% of the originally set intensity.
- modulation time (**TM** - possible choice: 2 sec. up to 60 sec)

Programming of Dense-Disperse-Stimulation according to Prof. Han (**CONT + BURST**) is not possible.

If you have chosen **“WORK/REST”**- rhythm you can set the following parameters:

- frequency (**Hz**)
- pulse width (**µS**)
- work time (**WORK**)
- rest time (**REST** = recovery time)
- increase / decrease ramp up time of pulse (**RU**)
- increase / decrease ramp down time of pulse (**Rd**)
- synchronous (**SY**) or alternating (**AL**) output of the channels A and B. If you have selected **“SY”** you can program **delay-time (DL)** from 0 sec. up to 5 sec. between the channels A and B.

6. By pressing the „S”button, you can set the next phase of the program. The active phase is indicated by a flashing arrow at the top of the screen.

The total number of arrows indicates the number of configured phases in this program.

You can create up to 5 phases of your own in the configurable program. If you want to reduce a configured program of 5 phases to less phases, you simply set the time of the first unwanted phase to zero.

7. The programming is completed by pressing the „P” button and the program is stored.

### Selectable frequency and pulse width settings:

Continuous-Mode: 1 to 150 Hz, 50 to 350 µS

Burst-Mode: 35 to 150 Hz, 50 to 200 µS

Modulated-Mode: 1 to 150 Hz, 50 to 350 µS

Work/Rest-Mode: 1 to 150 Hz, 50 to 450 µS

If frequency is > 100 Hz, pulse width is limited to 350 µS.

### The settings can be made as follows.

- Frequency: 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 15 Hz, and then increasing in 5 Hz steps until max. 150 Hz
- Pulse width: 50 to 450 µS in steps of 10 microseconds
- Intensity: 50-100% in steps of 10%
- Modulation time: 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 15, 20, 30, 40, 50, 60 seconds
- Work time: 2-99 seconds in steps of 1 second
- Ramp up time: 0,1-9,9 seconds in steps of 1 second
- Ramp down time: 0,1-9,9 seconds in steps of 1 second
- Rest time: 2-99 seconds in steps of 1 second
- Delay: 0,0-5,0 seconds in steps of 0,1 second
- Session time: 1-599 minutes in steps of 1 minute

## Programs

PRG 01 = Muscle rehabilitation, muscle strengthening, increasing resistance to muscular fatigue  
- high exertion

PRG 02 = Muscle rehabilitation, muscle strengthening - less exertion

PRG 03 = Relaxation and detoxification of the muscle, acceleration of the removal of toxic  
metabolic waste products, prevention of muscle cramps

PRG 04 = Reduction and/or elimination of pain

PRG 05 can be configured in up to 5 different phases

<b>Program</b>	<b>01</b>	<b>02</b>	<b>03</b>	<b>04</b>	<b>05</b>
<b>Mode</b>	<b>W/R</b>	<b>W/R</b>	<b>MOD</b>	<b>HAN</b>	<b>selectable</b>
Frequency	50 Hz	35 Hz	8-12 Hz	100/2 Hz	1-150 Hz
Pulse width	300 $\mu$ S	450 $\mu$ S	300-200 $\mu$ S	200/350 $\mu$ S	50-450 $\mu$ S
Running Time	20 min.	20 min.	30 min.	30 min.	1-599 min.
Ramp up	1 sec.	1 sec.			0,1-9,9 sec.
Ramp down	1 sec.	1 sec.			0,1-9,9 sec.
Work Period	10 sec.	3 sec.			2-99 sec.
Rest Period	20 sec.	15 sec.			2-99 sec.

**Abbreviations:**


AdJ:	Adjustment of intensity
AL:	Alternate function of the channels A and B
BURST:	The burst mode delivers 2 blocks of 9 impulses each per second
CONT:	In the continuous mode a permanent unchanged current is used
dA:	Day
dEL:	Delete
dL:	Delay between channels A and B
HI:	Upper value (high)
Hz:	Hertz is the measurement unit for the frequency of impulses per second.
LO:	Lower value
mA:	Milli-amp is the measurement unit for the electrical current intensity.
MOD:	In the modulated mode the frequency and/or the pulse width and/or the intensity changes in a predefined rhythm automatically.
Rd:	Ramp down
RU:	Ramp up
SEC:	Seconds
SY:	Synchronous function of the channels A and B
TI:	Time
TM:	Time of modulation-interval
W/R:	Work/Rest mode
µS:	Microsecond is the measurement unit for the length of every single impulse.

## LOCK-button and Statistics


The **Rehalito** offers the ability to lock one program in place so that other programs or settings cannot be selected or changed. This allows the doctor or therapist to ensure that the patient uses only the program that has been selected for him. The therapist can read information about the compliance of the patient on the display at any time he wants. This way it is possible to evaluate, whether the device was used according to the instructions of the doctor (quantity, intensity, etc.)

For this purpose the LOCK-button in the battery-compartment has to be pressed with a small pin once (the metal-pins of the electrode-cables are well suited). The lock symbol appears on the display and the program cannot be changed by pressing the „P“ -button anymore. Pressing the LOCK-button a second time allows reading the following statistics for the treatment in the display:

If the program was not locked, it is necessary to press the LOCK button twice to read the memory. The statistics are always recorded regardless of whether the unit is locked or unlocked. The statistics should be reset after 60 days of use or if the device is used by another person.

To get more detailed information press the + button of the left channel (A). You will see the data of the first session of the first day.


## LOCK-button and Statistics

The next sessions for the selected day are displayed - one after the other - by pressing the „S”- button. The next day can be obtained by pressing the + button of channel A again. Hold A+ or A- for three seconds to return to the global statistics screen.

If you require information about time and date of one single session please press any of the right +/--buttons (channel B). With these buttons you can scroll between session data, date and time.

The data stored in memory is retained after viewing, so that therapy-control is possible anytime. A maximum of 60 days of data can be stored with 5 sessions per day. When 60 days of data are reached, consecutive days will overwrite the earliest days. The day number will continue increasing (61, 62, ... up to max. 180) and will be displayed when the unit is unlocked until the **Rehalito** is reset.

If you press the P button for three seconds, you return to the home screen of the selected program. The device will not be locked. Press the lock button once to lock the device if it is required.

If you do not press „P” for 3 seconds, you can turn off the unit to leave it in lock mode when it is restarted.

The statistics should be reset before a new patient uses the device. Press LOCK to enter the statistics display and press „S” twice for 5 seconds each time. All settings in program PC05 will be reset as well. The **Rehalito** is set to its default state and all statistics are deleted.

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To achieve the best possible effect, one of the electrodes should be put on the motor point of the respective muscle if possible. Try to find out this point by slightly moving the electrode around. The other electrode is fastened on the remote muscle end. However, you also obtain a good effect if you put the electrodes respectively on the opposite muscle ends.

■ = red cable connection and ■ = black cable connection

**While stimulating, please keep this device in reach, to be able to finish the treatment at once, if current becomes uncomfortable.**

**Please ensure you follow the warnings and electrode placement advices on pages 23-25 of this manual.**

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In most cases it is very easy to find the suitable points for the placing of the electrodes. It is recommended to seek advice on electrode placement from your doctor or physiotherapist.

We suggest placing the electrodes on either side of the pain if no advice is available.

The distance between the electrodes should not be more then 20 cm. The polarity of the electrodes has no significance in TENS therapy.

Examples of possible placement-points are on the following pages.

The intensity should be set to a level that can be felt and then reduced slightly. The intensity setting should not cause visible muscle contraction.

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**Please ensure you follow the warnings and electrode placement advice on pages 23-25 of this manual.**

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









In accordance with European regulations [1], used electrical and electronic equipment may not be discarded in house hold or unsorted municipal waste. Use separate collection facilities. Contact your local government or council for information regarding the collection systems available. If electrical appliances are disposed of in landfills or dumps, hazardous substances can leak into the ground water and get into the food chain damaging the health and well being of people and animals.

[1] European Community directive 2002/96/EC on waste electrical and electronic equipment.

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
The battery should be replaced when it is drained or stops functioning. The battery supplied with the Rehalito is recyclable. Remove the battery from the unit and follow your local recycling guidelines.

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## Electromagnetic Compatibility

<b>Guidance and manufacturer's declaration – electromagnetic emissions</b>		
The Rehalito unit is intended for use in the electromagnetic environment specified below. The customer or the user of the Rehalito should ensure that it is used in such an environment.		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>
RF emissions CISPR 11	Group 1	The Rehalito uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Rehalito is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not Applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable	

<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
The Rehalito unit is intended for use in the electromagnetic environment specified below. The customer or the user of the Rehalito should ensure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment – guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	+/-6 kV contact +/-8 kV air	+/-6 kV contact +/-8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

## Electromagnetic Compatibility

<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
The <b>Rehalito</b> unit is intended for use in the electromagnetic environment specified below. The customer or the user of the <b>Rehalito</b> should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 V<sub>rms</sub> 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2,5 GHz</p>	<p>3 V</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the <b>Rehalito</b> unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p> <p><math>d=1.2\sqrt{P}</math>  <math>d=1.2\sqrt{P}</math> 80MHz to 800MHz  <math>d=2.3 \sqrt{P}</math> 800MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: right;">  </div>
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<p>A) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Rehalito unit is used exceeds the applicable RF compliance level above, the Rehalito unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Rehalito unit.</p> <p>B) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.</p>			

## Electromagnetic Compatibility

<b>Recommended separation distances between portable and mobile RF communications equipment and the Rehalito unit.</b>			
The <b>Rehalito</b> unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the <b>Rehalito</b> unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Rehalito unit as recommended below, according to the maximum output power of the communications equipment.			
<b>Rated maximum output power of transmitter W</b>	<b>Separation distance according to frequency of transmitter m</b>		
	<b>150 kHz to 80 MHz</b>	<b>80 MHz to 800 MHz</b>	<b>800 MHz to 2,5 GHz</b>
	$d = \left[\frac{3,5}{V_1}\right]\sqrt{P}$	$d = \left[\frac{3,5}{E_1}\right]\sqrt{P}$	$d = \left[\frac{7}{E_1}\right]\sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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Reserved for changes.

Printing or duplication of this manual - also in an abridged version - is forbidden and can only be allowed with a written approval of “MTR+ Vertriebs GmbH”.

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